SMALL CAPSULE. BIGIMPACT.

Celebrating 20 years of PillCam[™] capsule endoscopy





REVIEW ARTICLE



European guideline on indications, performance and clinical impact of ¹³C-breath tests in adult and pediatric patients: An EAGEN, ESNM, and ESPGHAN consensus, supported by EPC

Jutta Keller¹ | Heinz F. Hammer² | Paul R. Afolabi³ | Marc Benninga⁴ | Osvaldo Borrelli⁵ | Enrique Dominguez-Munoz⁶ | Dan Dumitrascu⁷ | Oliver Goetze⁸ | Stephan L. Haas⁹ | Bruno Hauser¹⁰ | Daniel Pohl¹¹ | Silvia Salvatore¹² | Marc Sonyi^{2,13} | Nikhil Thapar^{5,14} | Kristin Verbeke¹⁵ | Mark R. Fox^{11,16} on behalf of the European ¹³C-breath test group

Members of European 13C-Breath Test Group: Afolabi Paul; Altorjay, Istvan; Barbara, Giovanni; Basilisco, Guido; Baumann-Durchschein, Franziska; Belei, Oana; Benninga, Marc; Borrelli, Osvaldo: Churchev, Stanislav: Dominguez-Munoz, Enrique: Dumitrascu, Dan: Effenberger, Maria: Fox, Mark: Fürst, Stefan: Gasbarrini, Antonio: Goetze, Oliver: Haas, Stephan: Hammer, Heinz: Hammer, Johann; Hammer, Karin; Hauser, Bruno; Herszenyi, Laszlo; Keller, Jutta; Krznaric, Zeljko; Leja, Marcis; Lopetuso, Loris; Mion, Francois; Mulak, Agata; Nakov, Radislav; Nakov, Ventsislav; Pohl, Daniel; Reinisch, Sieglinde; Salvatore, Silvia; Shvet, Oleg; Simren, Magnus; Surdea-Blaga, Teodora; Tepes, Bojan; Thapar, Nikhil; Törnblom, Hans; Tutuian Radu; Verbeke, Kristin; Vogelsang, Harald; Vranesic-Bender, Darija; Wilder-Smith, Clive

This is an open access article under the terms of the Creative Commons Attribution-NonCommercial-NoDerivs License, which permits use and distribution in any medium, provided the original work is properly cited, the use is non-commercial and no modifications or adaptations are made.

© 2021 The Authors. United European Gastroenterology Journal published by Wiley Periodicals LLC on behalf of United European Gastroenterology.

¹Department of Internal Medicine, Israelitic Hospital, Academic Hospital University of Hamburg, Hamburg, Germany

²Department of Internal Medicine, Division of Gastroenterology and Hepatology, Medical University of Graz, Graz, Austria

³NIHR Southampton Biomedical Research Centre, University Hospital Southampton NHS Foundation Trust and University of Southampton, Southampton, UK

⁴Department of Pediatric Gastroenterology, Hepatology and Nutrition, Emma Children's Hospital, Amsterdam UMC, University of Amsterdam, Amsterdam, The

⁵UCL Great Ormond Street Institute of Child Health and Department of Gastroenterology, Neurogastroenterology and Motility, Great Ormond Street Hospital, London, UK

⁶Department of Gastroenterology and Hepatology, University Hospital of Santiago de Compostela, Santiago, Spain

⁷Clinica Medicala 2, Cluj-Napoca, Romania

⁸Department of Medicine II, Division of Hepatology, University Hospital Würzburg, Würzburg, Germany

⁹Department of Upper GI Diseases, Karolinska University Hospital, Stockholm, Sweden

¹⁰Department of Paediatric Gastroenterology, Hepatology and Nutrition, KidZ Health Castle UZ Brussels, Brussels, Belgium

¹¹Division of Gastroenterology and Hepatology, University Hospital Zürich, Zürich, Switzerland

¹²Pediatric Department, Hospital "F. Del Ponte", University of Insubria, Varese, Italy

¹³Clinic for General Medicine, Gastroenterology, and Infectious Diseases, Augustinerinnen Hospital, Cologne, Germany

¹⁴Department of Gastroenterology, Hepatology and Liver Transplantation, Queensland Children's Hospital, Brisbane, Australia

¹⁵Translational Research Center for Gastrointestinal Disorders, KU Leuven, Leuven, Belgium

¹⁶Digestive Function: Basel, Laboratory and Clinic for Motility Disorders and Functional Gastrointestinal Diseases, Centre for Integrative Gastroenterology, Klinik Arlesheim, Arlesheim, Switzerland

Correspondence

Jutta Keller, Israelitic Hospital, Orchideenstieg 14, 22297 Hamburg, Germany.

Email: j.keller@ik-h.de

Funding information UEG activity grant

Abstract

Introduction: ¹³C-breath tests are valuable, noninvasive diagnostic tests that can be widely applied for the assessment of gastroenterological symptoms and diseases. Currently, the potential of these tests is compromised by a lack of standardization regarding performance and interpretation among expert centers.

Methods: This consensus-based clinical practice guideline defines the clinical indications, performance, and interpretation of 13 C-breath tests in adult and pediatric patients. A balance between scientific evidence and clinical experience was achieved by a Delphi consensus that involved 43 experts from 18 European countries. Consensus on individual statements and recommendations was established if $\geq 80\%$ of reviewers agreed and <10% disagreed.

Results: The guideline gives an overview over general methodology of ¹³C-breath testing and provides recommendations for the use of ¹³C-breath tests to diagnose *Helicobacter pylori* infection, measure gastric emptying time, and monitor pancreatic exocrine and liver function in adult and pediatric patients. Other potential applications of ¹³C-breath testing are summarized briefly. The recommendations specifically detail when and how individual ¹³C-breath tests should be performed including examples for well-established test protocols, patient preparation, and reporting of test results.

Conclusion: This clinical practice guideline should improve pan-European harmonization of diagnostic approaches to symptoms and disorders, which are very common in specialist and primary care gastroenterology practice, both in adult and pediatric patients. In addition, this guideline identifies areas of future clinical research involving the use of ¹³C-breath tests.

KEYWORDS

breathtest, diagnosis, gastroenterology, gastroparesis, helicobacter pylori, liver cirrhosis, motility, pancreatic exocrine insufficiency, pancreatitis

INTRODUCTION

Breath tests are valuable, noninvasive diagnostic tests that are widely applied for the assessment of gastroenterological symptoms and diseases. ¹³C-breath tests provide the opportunity to diagnose *Helicobacter pylori (H. pylori)* infection, document gastric emptying time, monitor pancreatic exocrine and liver function, and have several additional potential gastroenterological and nongastroenterological applications.

Currently, the potential of breath testing is compromised by a lack of standardization regarding performance and interpretation among expert centers. This is highly relevant because modifications of the volume and/or composition of the test meal, of test performance and of the evaluation of data may markedly influence test results, diagnosis and thus, clinical usefulness of the investigation.

This consensus-based clinical practice guideline is needed within the gastrointestinal (GI) community to enhance pan-European harmonization of diagnostic approaches to symptoms and disorders, which are very common in specialist and primary care gastroenterology practice, both in adult and in pediatric patients. The guideline can add significantly to quality of investigation and, thus, the welfare of gastroenterological patients because it will allow a more rational approach to diagnostic evaluation and treatment. The guideline also aims to minimize disparities between health care systems across Europe, to facilitate cooperation between expert groups and the performance of multicenter clinical trials.

METHODS

The structured procedure, which was developed for the creation of this consensus-based clinical practice guideline, has previously been published. Briefly, this procedure was initiated by three representatives of the contributing societies (heads of guideline, JK, HH, MF) and started with formation of a representative core group of experts nominated from all participating societies and associations. This core group developed statements and recommendations, which were then submitted to reviewers in a three-stage Delphi voting process. The

heads of guideline and the core group members are listed as authors; the reviewers are listed as members of the European ¹³C-breath test group.

The following key questions were addressed in the guideline:

- 1. What is the role of ¹³C breath tests in the detection of *H. pylori* infection, and in the measurement of gastric emptying, pancreatic exocrine and liver function?
- 2. What are the general technical requirements and operating procedures for performance of ¹³C-breath tests, including preparation, dosage, breath sampling, technical analysis?
- 3. What are the reporting requirements?
- 4. Are there areas of disagreement and research priorities?

A systematic literature search with the appropriate key words using Medline/Pubmed and the Cochrane database was performed. We limited our search to studies performed in humans, which were published between 01 January 2000 and 25 July 2019. The resulting 446 references were assessed and allocated to the following topics: general methodology, ¹³C-urea breath tests (¹³C-UBT), ¹³C-gastric emptying breath tests (¹³C-GEBT), ¹³C-pancreatic function breath tests (13C-PFBT), 13C-liver function breath tests (13C-LFBT), and other 13C-breath tests. Statements and recommendations were developed based on these results, relevant consensus documents including those of participating societies (published after the year 2000),² and on pertinent literature known to members of the core group. Statements reflect key aspects and definitions but give no direct instructions on how to act, whereas recommendations advise when and how to perform individual ¹³Cbreath tests and how to report on breath test results. The wording used to indicate the strength of recommendation is detailed in Table 1.1

Of note, we aimed to develop recommendations on the specific indications for ¹³C-breath tests, whereas general indications for *H. pylori* testing, gastric emptying testing and monitoring of pancreatic exocrine or liver function were cited from current national and international guidelines.

Statements and recommendations were distributed among all reviewers via email for voting and commenting. All reviewers voted on all statements and recommendations according to the 6-point Likert scale given in Table 2^1 and gave comments in case of disagreement. "Agreement" was established if \geq 80% of reviewers voted A+/A AND < 10% D+/D.

Statements and recommendations, which did not receive agreement during the first Delphi round, were modified according to comments as previously described, and all modified statements and recommendations underwent further rounds of the Delphi process. Statements that achieved "agreement" after three Delphi rounds were accepted. Percentages for agreement and disagreement are given for all recommendations and statements. Recommendations additionally include quality of evidence (Q) and strength of recommendation (S). 1

The guideline is organized in a way that information, which will enable the reader to perform clinical tests, is summarized in statements, recommendations and respective tables. Background information detailed in the comments explains the choice of test protocols and further methodological aspects.

GENERAL 13C-BREATH TEST METHODOLOGY

Since their introduction into clinical practice in the 1970s,³ stable isotope ¹³C-breath tests have gained considerable importance and

TABLE 2 Six-point Likert-scale¹

Point	Description
A+	Agree strongly
Α	Agree with minor reservation
A-	Agree with major reservation
D-	Disagree with major reservation
D	Disagree with minor reservation
D+	Disagree strongly

TABLE 1 Descriptors of grading¹

Descriptor	Meaning	Wording
A-Strength high	Evidence or general accord that the procedure or statement is useful or effective. Further	has to be
	research is very unlikely to change our confidence in the estimate of effect	is to be
		shall
B—Strength moderate	Conflicting evidence or discordant opinions that the procedure or statement is useful or effective. The weight of evidence/opinion is in favor of utility. Further research is likely to have an important impact on our confidence in the estimate of effect and may change the estimate	should
C—Strength low	Conflicting evidence or discordant opinions that the procedure or statement is useful or effective. Further research is very likely to have an important impact on our confidence in the estimate of effect and is likely to change the estimate	could
D-Strength very low	Any estimate of effect is very uncertain	may

have been recommended for distinct diagnostic purposes by national and international guidelines and expert consensus papers. 4-7 This is due to the fact that they are not only reliable, but also noninvasive, relatively simple and safe diagnostic tools. ¹³C is a stable, nonradioactive carbon isotope with a natural abundance of about 1% of carbon isotopes. The use of non-radioactive stable-isotope tracers in biomedical experiments and diagnosis is generally considered ethically acceptable in humans at all ages. Toxicity of ¹³C has been examined in animals given amounts far in excess of those employed in the clinic. Up to sixty percent enrichment with ¹³C was achieved over prolonged periods of time without negative effects in adult animals and without signs of teratogenicity or embryotoxicity.8-10 Thus, the small requirements of ¹³C as a tracer in most clinical studies, particularly relative to its naturally high abundance, precludes any discernible risk of toxicity. Accordingly, 13C-enrichment of marker substances does not affect their tolerability. Depending on pharmacological properties of the marker substance itself, most clinically established ¹³C-breath tests can also be performed (repeatedly) in young children and pregnant women. This chapter gives important background information on methodological aspects and delineates general recommendations on test performance.

Statement 1.1 ¹³C-breath tests are used for investigation of a variety of gastrointestinal and liver functions and for diagnosis of Helicobacter pylori infection (100%, 0%).

In addition to the clinically established indications extensively discussed in this consensus report, several other tests have been developed for gastroenterological and other purposes (compare "Other tests").

Statement 1.2 The general principle requires that the digestive/metabolic process under investigation represents the rate-limiting step in the sequence of events leading to occurrence of ¹³CO₂ in the exhaled air (100%, 0%).

For instance, within the time frame used for testing, the presence of bacterial urease in the stomach determines whether 13 C-urea is metabolized leading to a selective increase in 13 CO₂-exhalation in *H. pylori* positive patients.

Statement 1.3 Most ¹³C-breath tests require sample collection and measurement over several hours. For selected indications (e. g., ¹³C-urea breath test in adults) measurements at two time points can be sufficient (100%, 0%).

Analysis of ¹³C-breath tests is frequently based on kinetic data of the ¹³CO₂ exhalation characteristics (e.g., gastric emptying tests) or on quantitative analysis of the whole metabolic process, that is, by analysis of cumulative ¹³CO₂ exhalation (e.g., pancreatic or liver function tests). Accordingly, several breath samples have to be collected at predefined intervals. Depending on the metabolic and

whole body distribution pathways of each ^{13}C -labeled substrate and the process under investigation, breath sampling over several hours may be necessary. For instance, protocols for pancreatic function testing usually require 4 to 9 h of sample collection at 15 or 30 min intervals. $^{11-14}$ On the other hand, ^{13}C -UBT in adults are usually based on two point measurements. 15,16 However, even for *H. pylori* testing, some experts calculate cumulative ^{13}C -exhalation over 1 h based on breath samples collected at 15 min intervals. The major advantage is that the patient's anthropometry is taken into account: a small and light person has a much higher measurement of δ over baseline value (DOB) than a very tall and heavy person for the same cumulative percentage of administered dose. In this way, the test can also be applied in children.

Statement 1.4 Isotope-ratio mass spectrometry is the reference method for measurement of the $^{13}CO_2$ concentration in the exhaled air (100%, 0%).

Statement 1.5 Isotope-selective nondispersive infrared spectroscopy can be used alternatively (91%, 0%).

At present, mass-spectrometry is the most accurate and efficient method for measuring carbon isotope ratios in exhaled breath (IRMS), but its application is restricted by the high cost of the equipment and operational complexity. Nondispersive infrared spectrometry (NDIRS) is the most widely used alternative method. Apart from lower costs, the devices are smaller, easier to handle and can be used on site, for example, in outpatient facilities. On the other hand, NDIRS measurements usually require higher sample volumes. Samples are frequently collected in aluminum bags (200–1300 ml) instead of 10 ml glass tubes as used for IRMS, which limits its use in large laboratories to which samples are delivered from distant sites.

Studies comparing the results of IRMS and NDIRS measurements of identical samples have shown comparable results for both methods. Most of these studies have investigated samples from ¹³C-UBT.¹⁷⁻²¹ Data show correlation coefficients of up to 0.999 for both analytical methods. A meta-analysis assessing the diagnostic accuracy of ¹³C-UBT in adult patients with dyspepsia showed no significant difference for studies reporting NDIRS or IRMS results.²²

Similarly, results of ^{13}C -octanoic acid breath tests (^{13}C -OABT), which measure gastric emptying, have been compared using IRMS and NDIRS. As expected, precision and repeatability of ^{13}C -measurements with NDIRS were inferior to IRMS. However, correlation coefficients for ^{13}C -exhalation and all gastric emptying parameters as computed on the basis of IRMS and NDIRS measurements were >0.98. Mean gastric emptying half time calculated using nonlinear regression (NLR) analysis was almost identical (87 \pm 39 min vs. 90 \pm 39 min). For the ^{13}C -methacetin breath test, a dynamic liver function test, molecular correlations spectroscopy, a method similar to NDIRS, showed comparable results to IRMS in Bland-Altman and correlation analysis. 24 Breath samples were

collected continuously via a nasal cannula with one sample analyzed about every 3 min.

Theoretically, these results should be transferable to other 13 C-breath tests since 13 CO $_2$ is the ultimate metabolic product, which is analyzed in all of these tests. However, caution may be necessary when applying tests, that only result in small increases in 13 C-exhalation, for example, 13 C-pancreatic function tests, for which highly accurate measuring devices are required. 25

Statement 1.6 The δ -value (‰) is the measuring parameter and is defined as the $^{13}\text{CO}_2/^{12}\text{CO}_2$ ratio in a given sample in comparison to the $^{13}\text{CO}_2/^{12}\text{CO}_2$ ratio in a reference material (97%, 0%).

Statement 1.7 Differences between δ -values obtained after application of a marker substance and the baseline δ -value (δ over baseline, DOB) are used for calculation of outcome parameters (97%, 0%).

Conventionally, 13 C-content of breath samples is expressed as δ % PDB units, that is, relative to the international standard, which originally was the calcium carbonate of the fossil Belemnitella of the Pee Dee formation (PDB) in South Carolina, USA. Zero δ ‰ corresponds to 1.112372% ¹³C atoms within CaCO₃. Thus, if a sample of carbon dioxide has a 13C/12C ratio which is less than that of the standard by 5 per mil, it is said to have a δ -value of -5%. ^{26,27} DOB values are the differences between δ-values obtained before (baseline) and after application of the marker substance. They reflect the increase in 13C-exhalation, which is the basis for calculation of test parameters. For most established ¹³C-UBT-protocols, the DOB value at 30 min is the relevant outcome parameter (compare below). For most other tests, based on DOB values and the assumption of a stable CO₂-production rate of 300 mmol per square meter of body surface per hour, the quantity of ¹³C appearing in breath per unit time is calculated. These data are usually expressed as percentage of ¹³C-dose administered.

Statement 1.8 13 CO $_2$ concentrations in samples collected using breath tubes for isotope-ratio mass spectrometry (IRMS) remain stable for at least 4 weeks so that measurement of breath samples can be delayed by this period (84%, 0%).

Statement 1.9 13 CO $_2$ concentrations in samples collected using aluminum bags for nondispersive infrared spectroscopy (NDIRS) remain stable for at least 72 h so that measurement of breath samples can be delayed by this period (93%, 0%).

For IRMS, stability of samples for a minimum period of 4 weeks has been demonstrated. A brief report even suggests that with 10 ml samples stored in glass tubes at room temperature in the absence of light, ¹³C-concentrations are stable for 8 months.²⁸ For NDIRS there are hardly any published data on sample stability. Mana et al. showed that samples are stable for 72 h.¹⁸ Personal

experience of the authors suggests longer stability (1-2 weeks by additional sealing of the rubber tubes of aluminum bags with gastight tapes).

Statement 1.10 Digestive and metabolic processes can be influenced by several factors including demographic parameters, fasting or fed state, composition and size of a test meal, physical activity, pre-existing diseases and drug intake (100%, 0%).

¹³C-exhalation from marker substances depends not only on the process under investigation (e.g., gastric emptying), but also on absorption of the marker substance and/or its metabolites, further (mostly hepatic) metabolism leading to production of ¹³CO₂, its transport to the lung and pulmonary excretion and distribution to other body compartments leading to a relevant loss of label, for example, into muscles or bone.²⁹ All of these functions as well as the process under investigation can be influenced by demographic and physiological parameters, concomitant diseases and drug intake.³⁰⁻⁴⁰ A study investigating the influence of clinical parameters on the results of ¹³C-OABT in more than 1200 patients has shown that ¹³C-exhalation was increased in women and correlated with age. Diabetes mellitus and inflammatory bowel disease were associated with decreased, and bacterial overgrowth and malignant disease with increased ¹³C-exhalation.³⁹

Figure 1 summarizes the general principle of ¹³C-breath testing.

What are the general prerequisites for performance of ¹³C-breath tests for clinical or research purposes?

Recommendation 1.1 ¹³C-breath tests performed for clinical reasons have to adhere strictly to standardized study protocols adequately validated in a representative patient population (100%, 0%; Q:C S: A).

Recommendation 1.2 For research projects, test parameters such as composition of the test meal or duration of breath sampling can be varied to evaluate the impact of the variation on test results (89%, 3%; Q:D S:B).

As discussed above, ¹³C-breath tests are indirect tests depending on several digestive and metabolic processes, which ultimately lead to exhalation of ¹³C-enriched breath. All intermediary steps as well as the process under investigation can be influenced by demographic and physiological parameters, concomitant diseases and drug intake. Moreover, alterations of marker substance, dose, or other components of the test meal/solution can markedly influence test results. Accordingly, it is of pivotal importance that tests are validated in representative patient populations and that clinicians adhere to strictly standardized study protocols. Research projects are required to better delineate the impact of variations on test results in order to optimize test procedures.




FIGURE 1 General principle of 13 C-breath testing. The 13 C-marked substrate is ingested, frequently with a specific test meal or test solution. Hepatic or bacterial enzymes release 13 CO₂, this requires prior intraluminal metabolization and/or intestinal absorption for most tests. 13 CO₂ is absorbed into the blood stream, transported to the lung, and exhaled. Breath samples are collected in glass tubes (mass spectroscopy) or aluminum bags (infrared spectroscopy) for measurement of 13 C/ 12 C-ratio. Delta over basal values are used for calculation of 13 C-excretion. * for specific liver function tests the 13 C-marked substrate is applied intravenously

Which dietary restrictions need to be observed before ¹³C-breath testing?

Recommendation 1.3 13 C-rich food, ingested before the test (e.g., corn, pineapple, broccoli, sugarcane) can increase the baseline δ-value and thereby compromise the measurements. Accordingly, they should be avoided at least 48 h before 13 C-breath testing (85%, 0%; Q:C S:B).

Isotopic fractionation—change in isotopic ratios between materials, due to the different rates at which various isotopes undergo chemical reactions—is a well-established phenomenon.²⁷ Carbon isotopes are strongly fractionated during photosynthesis, when plants metabolize carbon dioxide. Three types of photosynthesis occur in the plant world, commonly referred to as the C3, C4, and CAM pathways.²⁷ While most plants traditionally consumed in European diets perform C3 photosynthesis, leading to comparably lower ¹³C-content, other plants such as corn, pineapple, broccoli, and sugarcane are C4 plants with relatively higher ¹³C-abundance. Their consumption prior to a ¹³C-breath test increases basal ¹³C-exhalation, and further metabolization may influence ¹³C-exhalation over time. Avoidance for 48 h before breath testing is deemed satisfactory by most experts, while some recommend 72 h.

Which drugs and medical interventions need to be avoided before and during ¹³C-breath testing?

Recommendation 1.4 Drugs with potential influence on test results should be avoided before the test, unless essential long-term

medication is concerned or the effect of the drug on the digestive/metabolic process is to be determined (92% 0%; Q:D S:B).

Recommendation 1.5 Dialysis solutions and glucose infusions mostly contain glucose that originates from hydrolysis of maize starch naturally enriched in ¹³C and should therefore be avoided during ¹³C-breath testing (94%, 0%; Q:D S:B).

Drugs can influence test results by altering GI transit, absorption or (postabsorptive) metabolism of the ¹³C-labeled substrate. Accordingly, reliable performance of ¹³C-breath tests may require avoidance of specific drugs as discussed for the individual tests below. However, if long-term treatment is mandatory, drug avoidance is not always reasonable because it confounds the normal clinical situation of the patient. Moreover, ¹³C-breath tests can be used for monitoring drug effects. For instance, the ¹³C-mixed triglyceride breath test (¹³C-MTGBT) has been used to monitor improvement of lipid absorption with enzyme replacement therapy in pancreatic exocrine insufficiency (PEI).¹¹

Dialysis solutions and glucose infusions mostly contain glucose that originates from hydrolysis of maize starch, so that they are naturally enriched in ¹³C. Accordingly, they may confound test results and should also be avoided.³⁷

Is physical activity allowed during ¹³C-breath testing?

Recommendation 1.6 Physical activity alters gastrointestinal transit of orally administered substrates and markedly increases CO₂

production. Therefore, physical activity has to be avoided during ¹³C-breath testing (100%, 0%; Q:C S:A).

Even moderate physical activity such as walking roughly doubles energy expenditure compared with sedentary subjects and has corresponding effects on endogenous CO₂ production. ^{41,42} In addition, exercise leads to a shift toward oxidation of nonlipid components. This increases ¹³C-exhalation because the lipid molecules in the body contain substantially lower concentrations of ¹³C than the nonlipid molecules, due to fractionation processes during lipid synthesis. ³⁵ Moreover, ¹³C-breath tests using orally applied marker substances depend on GI transit which is accelerated by moderate exercise, while strenuous exercise has opposite effects. ⁴³ Accordingly, it has been shown that physical activity during ¹³C-OABT markedly alters ¹³CO₂-exhalation in healthy volunteers as well as respective normal values. ³⁶

In summary, physical activity has profound and complex effects on ¹³C-exhalation and breath test results. For standardization purposes under clinical conditions, patients must be asked to strictly avoid physical activity during tests.

Which concomitant diseases may influence ¹³C-breath test results?

Recommendation 1.7 Disturbances of gastrointestinal motor and secretory function, hepatic and pulmonary function can generally affect the time course and/or amount of ¹³C-exhalation. This has to be taken into account for performance and interpretation of ¹³C-breath tests (100%, 0%; Q:B S:A).

Concomitant diseases affecting digestive and metabolic processes which ultimately lead to exhalation of ¹³C-enriched breath may confound breath test results. In particular, major disturbances of GI transit and absorption, hepatic and lung function have to be considered. However, the influence of concomitant diseases on many ¹³C-breath tests appears to be small:

Even in patients with nonalcoholic steatohepatitis or liver cirrhosis (\sim 50% Child-Pugh score C), 13 C-octanoic acid metabolism was found to be normal $^{44-46}$ and the 13 C-OABT correlated well with scintigraphy in patients who were critically ill. 47

Disturbances of $^{13}\text{CO}_2$ -exhalation in patients with lung diseases are only expected in patients with very severe disease and reduced CO_2 -diffusion capacity. Data on potential impairment of ^{13}C -breath tests in such cases are not available in the literature.

Since digestion of dietary lipids by pancreatic lipase cannot occur before the meal has entered the duodenum, results of the ¹³C-MTGT are influenced by and correlate with gastric emptying parameters. ¹⁴ This predominantly applies to early postprandial ¹³C-exhalation rates (1–3 h), which have insufficient specificity for detection of PEI. However, this problem can be overcome by prolongation of the breath sampling period so that gastric emptying is of minor importance for ¹³C-exhalation.

Vice versa, ¹³CO₂-excretion following administration of ¹³C-octanoic acid was similar in patients with and without PEI, including those with severe exocrine insufficiency characterized by overt steatorrhea. ^{39,48,49}

¹³C-UREA BREATH TEST

H. pylori is a common bacterial pathogen responsible for substantial GI morbidity worldwide. In addition to causing inflammatory gastroduodenal alterations, *H. pylori* is the major risk factor for gastric cancer development and is associated with various other, partly non-GI diseases. Table 3 summarizes important indications for *H. pylori* testing in adults and children recommended by current European guidelines.^{4,5}

There are differences in the approach to *H. pylori* infection between adults and children. Thus, the European Society of Pediatric Gastroenterology, Hepatology, and Nutrition (ESPGHAN) explicitly recommends against a test and treat strategy in children, as well as against *H. pylori* testing as part of the initial investigation in children

TABLE 3 Indications for *H. pylori* testing in adults and children^{4,5}

High	Suspicion/evidence of peptic ulcer disease, atrophic gastritis, gastric adenocarcinoma, MALT (mucosa-associated lymphoid tissue) lymphoma
	Test-and-treat strategy for uninvestigated dyspepsia
	Exclusion of H. pylori gastritis before reliable diagnosis of functional dyspepsia
Moderate	Aspirin and NSAIDs users with a history of peptic ulcer
Low	Unexplained iron deficiency anemia, idiopathic thrombocytopenic purpura, vitamin B12 deficiency
High	Suspicion/evidence of peptic ulcer disease
	Monitoring of outcome of eradication therapy
Low	Chronic immune thrombocytopenic purpura
	Moderate Low High

with iron deficiency or as part of the investigations of causes of short stature. Moreover, these guidelines recommend against *H. pylori* testing in children with functional abdominal pain disorders.⁵

For the diagnosis of H. pylori infection, histology (\pm culture) from biopsy samples is regarded as the reference standard. Several other tests are available, including endoscopic biopsy followed by rapid urease testing and noninvasive methods like serology, fecal antigen tests and the 13 C-UBT. The principle of the 13 C-UBT relies upon the capacity of H. pylori, when present in the stomach, to hydrolyze orally administered 13 C-urea to produce 13 CO₂, which diffuses into the blood, is transported to the lungs and exhaled, so that it can be detected in breath samples. 13 C-urea is innocuous and can be safely administered repeatedly, including in children and pregnant women.

According to the Maastricht V/Florence Consensus report, ¹³C-UBT is regarded as the best approach to the diagnosis of *H. pylori* infection in the context of a "test-and-treat strategy" in adults because of its high sensitivity and specificity, and excellent performance.⁴ This is confirmed by a recent Cochrane review showing superior diagnostic accuracy compared with other non-invasive tests.⁵⁰

Stool antigen tests may be less acceptable in some societies but also have a high sensitivity and specificity, provided a monoclonal antibody-based ELISA is used.⁴ Some serology tests have high sensitivity and specificity, but these tests may perform differently in different geographic locations according to the antigenic composition of the circulating strains. Thus, the Maastricht V/Florence Consensus recommends that only locally validated tests should be used.⁴ Moreover, due to the slow decrease of serum antibody levels, they are inadequate for assessment of *H. pylori* eradication after treatment.

In clinical practice when there is an indication for endoscopy, and there is no contraindication for biopsy, the rapid urease test (RUT) is recommended as a first-line diagnostic test. In the case of a positive test, it allows immediate treatment.⁴

When should ¹³C-UBT be utilized?

Recommendation 2.1 The ¹³C-urea breath test is to be considered as a noninvasive alternative for all indications for Helicobacter pylori testing if endoscopy is not required or if biopsies are contraindicated (92%, 0%; Q:A S:A).

Recommendation 2.2 The ¹³C-urea breath test is a preferred option for confirmation of Helicobacter pylori eradication in adults and children. It has to be performed at least 4 weeks after completion of therapy (89%, 0%; Q:A S:A).

Meta-analyses confirm that 13 C-UBT achieves \geq 95% sensitivity and specificity in both adults and children. 50,52 Accordingly, it is a highly accurate test and generally appropriate for all indications for *H. pylori* testing. However, in clinical practice, when there is an

indication for endoscopy, and there is no contraindication for biopsy, the RUT is recommended as a first-line diagnostic test.⁴ Some national guidelines suggest two positive tests for reliable diagnosis, except in cases with a very high a priori likelihood of *H. pylori* infection, for example, duodenal ulcer,⁵³ and ¹³C-UBT can be readily performed. Moreover, ¹³C-UBT is regarded as the best option for confirmation of *H. pylori* eradication,^{4,50} with stool antigen tests being an alternative. Testing has to be delayed for at least 4 weeks after the end of therapy, otherwise it may lead to false negative results (compare patient preparation).

How should ¹³C-UBT be performed?

Recommendation 2.3 If commercially available kits are used for ¹³C-urea breath test, manufacturers' instructions regarding preparation of test solution and test performance have to be followed (97%, 0%; O:B S:A).

Recommendation 2.4 If the test is prepared on site, investigators have to adhere to well established test protocols (for examples compare Table 4) (97, 0%, 0%; Q:B S:A).

As discussed above, alterations of marker substance, dose or other components of the test meal/solution can markedly influence test results. Hence, it is of pivotal importance that tests are validated in representative patient populations and that clinicians adhere to strictly standardized study protocols. Since the original description of the ¹³C-UBT by Graham et al.,⁵⁷ several changes of the test protocol have been proposed affecting the dose of ¹³Curea, type of test meal/solution, time of breath collection, cut-off values, and measuring device.⁵⁸ Logan et al. first introduced a modified version with only one breath sample to be analyzed, however this was pooled from several samples collected over 30 min at 5 min intervals.⁵⁹ Currently, most studies report DOB values of a single sample collected 30 min after application of ¹³Curea with a threshold of >4% for diagnosis of H. pylori infection. At this threshold the summary sensitivity (95% confidence interval [CI]) and specificity (95% CI) from 10 studies (958 participants) were 0.95 (95% CI: 0.79-0.99) and 0.95 (95% CI: 0.87-0.98).⁵⁰ A minority of studies used sampling periods of 10 or 20 min and thresholds between DOB>3‰ and DOB>6‰.50 13C-urea was usually applied with a citric acid solution; however, orange juice⁵⁴ or semiliquid meals have also been used.55 Examples for wellestablished test protocols are given in Table 4.

How should patients prepare for the test?

Recommendation 2.5 Ideally, adult patients should have fasted overnight. If this is not feasible, a fasting period of 4 to 6 h is sufficient (97%, 0%; Q:C S:B).

TABLE 4 Established test protocol for ¹³C-UBT in adults and children

Reference		Test solution	Breath sampling	Cut off	Validity	Remarks
Leodolter 1999 ¹⁵	Adults	75 mg ¹³ C-urea dissolved in 200 ml 0.1 mol/L citric acid solution (~4 g/200 ml water) with two tablets artificial sweetener	Before and 30 min after ingestion	≥4‰	SENS 95% SPEC 98% ACC 97%	Equal performance as test with citric acid solution 10 min before marker ingestion; With 200 ml orange juice instead of citric acid lower SENS (88%) with equally high SPEC (100%) ⁵⁴ ; with semiliquid meals longer test duration required ⁵⁵
Elitsur et al. 2009 ⁵⁶	Children	75 mg ¹³ C-urea and 2 g citric acid dissolved in 4 ounces (~120 ml) potable water	Baseline and 15 min after ingestion	≥2.4‰	SENS 97.9% SPEC 96.1% PPV 90.4% NPV 99.2%	Test performed best in children aged >6 years; in ages 2-5 calculation of urea hydrolysis rate can lead to higher SENS and SPEC compared with DOB values

Abbreviations: ACC, accuracy; DOB, δ over baseline; NPV, negative predictive value; PPV, positive predictive value; SENS, sensitivity; SPEC, specificity; UBT, urea breath test.

Recommendation 2.6 Children of all ages should have fasted for 4 h (100%, 0%; Q:B S:B).

Theoretically, food in the stomach may dilute the marker substance and thereby impair the contact between the infected mucosa and the substrate leading to lower ¹³CO₂ production and decreased sensitivity of the test. Although this aspect is still controversial, some studies suggest that fasting before the ¹³C-UBT indeed improves accuracy of test results. ^{18,60} Accordingly, it appears to be prudent to perform the ¹³C-UBT in fasting conditions. ⁶⁰ In adults, fasting overnight appears reasonable. If this is not feasible, a fasting period of 4-6 h is deemed sufficient. In young children, long fasting periods may be problematic. However, sensitivity of the ¹³C-UBT was markedly reduced to about 50% when children were fed a meal immediately before the test. ⁶¹ Therefore, a fasting period of 4 h is suggested, including in young children.

Recommendation 2.7 Prior to the 13 C-urea breath test patients have to abstain from: proton pump inhibitor therapy for ≥ 2 weeks; antibiotic therapy (including Helicobacter pylori eradication) for ≥ 4 weeks (97%, 0%; Q:A/B S:A).

Recommendation 2.8 Antacids can be allowed before the ¹³C-urea breath (86%, 0%; Q:B S:B).

Proton pump inhibitors (PPI) need to be discontinued because they decrease the load of *H. pylori* leading to false-negative results on several tests including ¹³C-UBT. ^{4,62,63} A 7-day withdrawal has been shown to be sufficient in most patients; however, as a precaution, 14 days are recommended. ⁴ H₂-receptor antagonists slightly decrease sensitivity of ¹³C-UBT for up to 2 weeks. ⁶⁴ Topical antacids do not affect sensitivity ⁶² so that they can be allowed before testing. By contrast, antibiotics, including those used for eradication therapy, and bismuth compounds need to be discontinued for 4 weeks to allow an increase of a detectable bacterial load. ⁴

How should test results be reported?

Recommendation 2.9 To allow for reliable interpretation of test results, the following parameters should be reported: marker dose and test solution; test result including normal values and interpretation (Helicobacter pylori negative/positive) (95%, 0%; Q:D S:B).

Recommendation 2.10 The test report could be complemented by including clinical characteristics of the patient, last use of proton pump inhibitor and the exact test protocol including equipment used for breath sampling and analysis (84%, 0%; Q:D S:C).

Experts agree that a minimum of information on methodology (marker dose and test solution) and test results (e.g., DOB-value) including normal values and interpretation (*H. pylori* positive/negative) are required for medical personnel not involved in the testing to reliably interpret individual findings. Clinical characteristics of the patient including last use of PPI and further methodological information may further facilitate assessment of reliability of test results and choice of clinical consequences.

¹³C-GASTRIC EMPTYING BREATH TESTS

Gastric dysmotility can manifest as rapid gastric emptying with dumping syndrome (even in the absence of upper GI surgery) or delayed gastric emptying with symptoms of gastroparesis. The latter applies to the majority of affected patients and is typically associated with nausea, vomiting, early satiety, postprandial fullness, upper abdominal pain, and bloating in adults and children. 6,65–69 Anorexia and weight loss are further frequent symptoms. Children with gastroparesis experience more vomiting while adolescents with gastroparesis report more nausea and abdominal pain. 68,70 There is general consensus that the diagnosis of gastroparesis requires objective evidence of clearly delayed gastric emptying in symptomatic

patients.^{6,66} However, prior to gastric emptying testing, the exclusion of mucosal or structural disorders such as inflammatory or malignant diseases as the underlying cause of symptoms is required.⁶ Specific indications for gastric emptying testing in adults as suggested by international guidelines and expert consensus papers are given in Table 5.^{6,65,71} For pediatric patients, there are no generally accepted guideline recommendations on gastric emptying testing. The recommendations given in Table 5 are derived from a recent review of the literature, which shows that nausea, vomiting, and abdominal pain are the most common symptoms in children, while early satiety, postprandial fullness, bloating, and weight loss occur less frequently in pediatric gastroparesis.⁷²

When should ¹³C-GEBT be utilized?

Recommendation 3.1 ¹³C-gastric emptying breath tests are to be regarded as an established alternative to scintigraphy for measurement of gastric emptying velocity (92%, 0%; Q:A S:A).

Scintigraphy is the reference standard for measurement of gastric emptying. However, while there is a consensus report recommending a standardized protocol in adults in the United States, ⁷³ no European consensus exists on the type of test meal and duration of data acquisition. Likewise, no consensus exists for a standard gastric emptying scintigraphy in pediatrics. However, recent studies have provided confirmation that extending studies from 2 to 4 h increases the diagnostic yield and should be the standard in children and adolescents as it is in adults for measurement of solid gastric emptying. ⁷³⁻⁷⁵ For liquids, 2 h are probably sufficient and early gastric emptying has to be accounted for. ⁷⁶

Several protocols for 13 C-based gastric emptying tests have been successfully validated in comparison with scintigraphy. The medium-

chain fatty acid, ¹³C-octanoic acid,^{77,78} or the edible blue-green algae, ¹³C-*Spirulina platensis*⁷⁹ are typically used to label solids; ¹³C-acetate is used for liquids.^{80,81} On delivery to the duodenum, the ¹³C-containing substrate is either absorbed directly (octanoic acid, acetate) or digested and then absorbed (*Spirulina*). Subsequently, it is metabolized in the liver, and finally excreted by the lungs as ¹³CO₂.⁶

The first use of the ¹³C-OABT in adults was published by Ghoos et al. in 1993.⁷⁷ Gastric half emptying time (T½) was assessed by NLR analysis and corrected for the expected delay caused by postgastric processes (absorption, metabolization, exhalation). The authors observed an excellent correlation between ¹³C-OABT parameters and parameters obtained from simultaneous scintigraphy (R = 0.89 for T½ scintigraphy vs. breath test). Sensitivity and specificity, positive and negative predictive values of breath test parameters for delayed gastric emptying were ≥94%.⁷⁷ Results of the ¹³C-OABT also closely correlate with those of scintigraphy in other studies in adults and children. 78,82-84 Following this initial study, all kinds of variation of the test meal and the test protocol have been described, depending on cultural differences and practical considerations (e.g., muffins, pancakes, rolls). In addition, different mathematical analysis methods have been proposed^{77,80,85-88} (also compare Table 6). Pancakes marked with ¹³C-octanoic acid are an acceptable and palatable solid test meal for children but cannot be used in case of allergy to egg, milk or wheat, in coeliac disease or very young children. 78,88,89

The lack of standardization is of concern because it makes comparison of values between different laboratories difficult. On the other hand, it allows flexibility in research projects to measure gastric emptying of very different test meals and evaluate the impact of composition on gastric emptying.

Results of the 13 C-Spirulina-GEBT also show high concordance (R=0.86) with scintigraphic data. 90 The protocol is exactly defined and has been validated in a large group of healthy volunteers and

 TABLE 5
 Indications for gastric emptying testing in adults and children 6,7,65,71,72

Grade of recommendat	ion	
Adults	High/moderate	Symptoms suggestive of gastroparesis* without evidence of mucosal or structural disease explaining these symptoms (*nausea, vomiting, early satiety, postprandial fullness, bloating, upper abdominal pain)
		Unexplained impairment of blood glucose control in patients with diabetes mellitus, even in the absence of abdominal symptoms (because of the central role of gastric emptying for regulation of postprandial glycemia)
	Moderate	To support the diagnosis of dumping syndrome
	Low	Severe gastroesophageal reflux disease unresponsive to acid suppressants (particularly before fundoplication); systemic sclerosis; after lung transplantation; Parkinson disease; generalized GI motility disorders; patients under consideration for intestinal or colonic surgery or transplantation because of motility disorders
Children	High/moderate	Most common GI symptoms located in the upper GI tract suggestive for gastroparesis: Nausea, vomiting and abdominal pain
	Low	Less frequent GI symptoms located in the upper GI tract suggestive for gastroparesis: Early satiety, postprandial fullness and bloating

Abbreviation: GI, gastrointestinal.

TABLE 6 Established test protocols for ¹³C-GEBT in adults and children^a

Reference	Age	Estimated parameter Test meal	Test meal	Breath sampling	Endpoints and normal values	Validity	Remarks
Ghoos 1993 ⁷⁷	Adults	Solid GE	Two slices of white bread, 5 g butter, 200 ml water, omelet made from one egg, yolk doped with 91 µg (= 100 µl) ¹³C- octanoic acid	Samples at baseline (preferentially taken as duplicate), further samples at 15 min intervals up to 4 h pp	T ½ (mean ± 25D): 28- 116 min ^c	T %: SENS 95% SPEC 94% PPV 94% NPV 94%	Evaluated against SCINTI in HC ($N = 16$) and patients with dyspepsia ($N = 20$), normal values from 42 HC (NLR model), no test kit commercially available, other groups report slightly different T ½ normal values using same protocol: 50–150 min ⁹⁹
Szarka 2008 ⁹⁰	Adults	Solid GE	Freeze-dried scrambled eggs mix containing 100 mg ¹³ C-Spirulina platensis, six saltine crackers, and 180 ml of water	Breath samples at baseline, on completion of the meal and at 45, 90, 120, 150, 180, and 240 min pp	kPCD values at 45, 150, and 180 min provide strongest concordance with scintigraphy for accelerated and delayed GE	Delayed GE: SENS 89%, SPEC 80% Accelerated GE: SENS 93%, SPEC 80%	38 HC and 129 patients with clinically suspected delayed GE, normal T ½ according to SCINTI (10th-90th percentile): 52-86 min,FDA approved, CE marked, commercially available in the United States, only
Bertram 2014 ¹⁰⁰	Adults	Liquid GE	150 mg ¹³ C-acetate dissolved in 200 ml water with 10 g lacutose	Breath samples at baseline, at 5 min intervals for first hour, at 15 min intervals for second hour pp	Time of maximal ¹³ C-exhalation ^b : (P10-P90): 15-40 min	Time of maximal ¹³ C-exhalation and T ½ SCINTI: R = 0.88, p < 0.005 in validation study by Chew 2003 ⁸⁷	22 HC, lactulose used for simultaneous measurement of liquid gastric emptying and small bowel transit by H ₂ -breath test, time to maximal ¹³ C-exhalation in HC identical with validation study ³⁷ in 10 HC which used 15 g glucose instead of lactulose
Van Den Driessche 1999 ¹⁰¹	29 healthy premature and term infants gestational age 2741 weaks, postnatal age 7-74 days	Liquid GE	Group 1: 50 ml expressed breast milk, Group 2: 50 ml infant formula (Nutrilon premium®) (33.5 kcal), each with 50 µl ¹³ C-octanoic acid	Breath sample at baseline, further samples at 5 min intervals for 30 min, then at 10 min intervals up to 4 h pp	T ½ (mean, range): group 1 = 47, 16-86 min ^c ; Group 2 = 65, 27- 98 min ^c		
							(Continues)

TABLE 6 (Continued)

Reference	Age	Estimated parameter Test meal	Test meal	Breath sampling	Endpoints and normal values	Validity	Remarks
Hauser 2016 ¹⁰²	133 healthy children mean 9 years, range 1–17 years	Liquid GE	200 ml INZA® milk-drink Breath samples at (skimmed milk) baseline and at (112 kcal) with 50 mg intervals for 4C (body weight 10- then at 10 min 30 kg) or 100 mg intervals up to (>30 kg) ¹³ C-acetate	Breath samples at baseline and at 5 min intervals for 40 min, then at 10 min intervals up to 3 h pp	T ½ (mean ±2SD): 55– 109 min ^b Normal values: Percen- tiles according to age Delayed gastric emptying defined as T ½ > P90 Rapid gastric emptying defined as T ½ < P10	T ½ 13 C-ABT and T ½ SCINTI: R = 0.604 (ρ = 0.0006) Reproducibility tested in 21 healthy children: CV T ½ (median, range) = 8.3, 1.6%-16.2%	Comparison with scintigraphy in 21 children with upper GI symptoms
Eradi 2006 ¹⁰³	25 healthy children mean 7.8 years, SD 0.3 years, range 5-10 years	Solid GE	30 g chocolate crispy cake (147 kcal) with 100 mg ¹³ C-octanoic acid	Breath samples at baseline, further samples at 15 min intervals up to 4 h pp	T ½ (mean ± 2SD): 44- 155 min ^b	Comparison 13 C-OABT and gastric emptying scintigraphy:1 12 C-OBT and T 12 -SCINT!: $R=0.69~(p<0.01)$,
Hauser 2016 ⁷⁸	120 healthy children mean 9 years, SD 4 years, range 1-17 years	Solid GE	One pancake (17 g wheat Breath sample at flour, 7 g sugar, one egg samples at 15 yolk, 40 ml semiskimmed milk, 5 g margarine) + 5 g sugar + 100 ml water (230 kcal) with 50 µl	Breath sample at baseline further samples at 15 min intervals up to 4 h pp	T ½ (mean ± 2SD): 50–266 min ^b Normal values: Percentiles according to age: Delayed gastric emptying defined as T ½ > P90 Rapid gastric emptying defined as T ½ < P10	T ½-13C-OABT and T ½- SCINTI: R = 0.748 (p < 0.0001) Reproducibility tested in 19 healthy children: CV T ½ (median, range) = 13.3, 2.56%-29.6%	Comparison with scintigraphy in 19 dyspeptic children

Abbreviations: ¹³C-ABT, ¹³C-acetate breath test; ¹³C-OABT, ¹³C-octanoic acid breath test; CV, coefficient of variance; GE, gastric emptying; GEBT, gastric emptying breath test; GI, gastrointestinal; HC, healthy controls; RPCD, percent dose excreted × 1000; NLR, nonlinear regression model; NPV, negative predictive value; P10, 10th percentile; P90, 90th percentile; pp, postprandial; PPV, positive predictive value; SCINTI, scintigraphy; SENS, sensitivity; SPEC, specificity; T½, gastric half emptying time.

^aOnly studies with $N \ge 20$.

^bCalculated breath test data not corrected for scintigraphy.

^cScintigraphic equivalent values or breath test data corrected for scintigraphy according to Ghoos et al. 1993,⁷⁷

patients.⁹⁰ The test was approved by the FDA for evaluation of gastric emptying in 2015 (https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpma/pma.cfm?id=P110015), and it is also CE marked according to the producer's information (https://sercon-instruments.com/wp-content/uploads/2017/06/019-gebt.pdf). However, it has not been marketed in Europe, so far.

¹³C-acetate has been used as a marker for acaloric and caloric liquids such as formula diets. ^{81,87,91–93} In young children or children dependent on gastrostomy feeding, the test meal often consists of milk or milk-based formulas. ⁹⁴ Apple juice has also been used. ³² Results of the test closely correlate with those of scintigraphy in adults and children. ^{81,87,91–93} Delayed gastric emptying of solids usually precedes disturbances in gastric emptying of liquids. ⁹⁵ Therefore, tests of solid gastric emptying are supposed to have a higher sensitivity. However, gastric emptying of liquids can be abnormal in patients with normal gastric emptying of solids. ⁹⁶ Moreover, liquid test meals are obviously more useful in young children ^{81,94} and probably also for confirmation of rapid gastric emptying in patients with suspected dumping syndrome. ^{6,97}

Intra-individual and inter-individual variabilities of all ¹³C-GEBTs are high in adults and children, but similar to variations observed with scintigraphy^{85,88-90,94,98} and, therefore, reflect day-to-day physiological variability in gastric emptying.⁶

How should ¹³C-GEBT be performed?

Recommendation 3.2 ¹³C-GEBT (gastric emptying breath test) performed for clinical reasons have to adhere strictly to standardized study protocols adequately validated in a representative patient population. This refers to preparation of the test meal as well as test performance and evaluation of test results (for examples compare Table 6) (97%, 0%; Q:C S:A).

As discussed above, ¹³C-GEBT are indirect tests that involve multiple steps and are prone to influences caused by demographic, physiological, and other parameters.^{39,99,104} It has been hypothesized that ¹³C-GEBT might be inaccurate in conditions associated with substantial malabsorption, liver, or lung diseases, though this is not substantiated by clinical studies.^{39,44,47} However, it is still important that tests are validated not only in healthy adults and children, but also in the relevant patient population.^{77–79,81,84}

Moreover, alterations of marker substance, dose or other components of the test meal markedly influence test results. For instance, a larger labelled test meal will result in higher normal ranges for $T_2^{1/2}$ and gastric lag time. For $^{1/2}$ C-acetate an interaction has been demonstrated between the rate of $^{1/2}$ C-delivery to the duodenum and $^{1/2}$ C-recovery in breath.

Different mathematical models have been developed for analysis of gastric emptying curves derived from breath tests, in particular the NLR model, the generalized linear regression (GLR) model and the Wagner–Nelson method.^{77,79,106} Cumulative ¹³CO₂-excretion over

time is inversely related but analogous to the scintigraphic gastric emptying curve. However, ¹³CO₂-excretion does not only depend on gastric emptying velocity but also on postgastric absorption and metabolisation of the substrate and CO₂-exhalation rates. For this reason, Ghoos et al. developed the original NLR model.⁷⁷ According to this model, T½ indicates the time at which half of the ¹³CO₂ is excreted, relative to the cumulative excretion when time is infinite. Accordingly, results are determined by the shape of the exhalation curve, independent of absolute ¹³CO₂-excretion. Measurements are usually performed over 4 h with breath samples at 15 min intervals.

The GLR model published by Lee et al. proposed a minimum number (N=3) of breath samples at pre-specified times during the 3 h postprandial period to mathematically predict the gastric emptying endpoints measured by simultaneous scintigraphy. Results reflect absolute $^{13}\text{CO}_2$ -excretion. A similar model with breath samples obtained upon completion of the meal and then at 45, 90, 120, 150, 180, and 240 min postprandially was suggested by the same group and is used for analysis of the test commercially available in the USA.

The Wagner–Nelson method has been suggested for analysis of $^{13}\text{C-GEBT}$ by Sanaka et al. 106 This method has been developed to describe the entrance of ingested drugs into the venous system based on its urinary excretion data. When applied on breath tests, it describes the manner in which $^{13}\text{CO}_2$ appears in the venous system based on pulmonary $^{13}\text{CO}_2$ excretion data. 106,107 It is used less frequently, and, similar to the NLR model, the $^{13}\text{CO}_2$ -exhalation curve must exhibit the decreasing portion during the sampling period for correct estimation of gastric emptying parameters. 108 Accordingly, breath sampling has to be routinely performed for 4 h and potentially longer in gastroparetic patients.

Given that in Europe there is no standardized, well-validated test kit commercially available, and that tests are usually prepared on site, it is of pivotal importance that clinicians adhere strictly to standardized study protocols including established analysis methods. Especially in children, a large variety of test meals have been explored in accordance with the variable requirements of different age groups. $^{78,101-103,109-114}$ However, several of these studies were performed in small patient groups. Examples of well validated 13 C-GEBT protocols in adults and children (studies with $N \geq 20$) are given in Table 6.

How should patients prepare for the test?

Recommendation 3.3 Before and during the test, precautions as described in General Methodology (avoidance of ¹³C-rich food, avoidance of physical activity, and ¹³C-rich infusions during test) have to be observed (100%, 0%; Q:C S:A).

Recommendation 3.4 Drugs with potential influence on gastrointestinal transit should be avoided before the test, unless essential long term medication is concerned or the test is performed to monitor the drug effect on gastric emptying (100%, 0%; Q:C S:B).

Recommendation 3.5 Adult patients, older children, and adolescents have to be fasted overnight (94%, 0%: Q:C S:A).

Recommendation 3.6 A shorter fasting period can be sufficient in very young children (92%, 0%; Q:C S:B).

Dietary and other restrictions, which generally apply before and during ¹³C-breath testing, have been explained in General Methodology. Physical activity, in particular, has to be avoided during ¹³C-GEBT, not only to standardize CO₂-production but also because physical activity influences gastric emptying velocity. ^{43,115} Tests should preferentially be performed in the sitting position since the supine position may be associated with slower gastric emptying. ¹¹⁶ Experts agree that solid test meals should be consumed within 10–15 min, liquids within 5–10 min.

Adults, adolescents and older children are required to fast overnight prior to breath testing^{77,79,88}; more precisely, a fasting period of ≥12 h is recommended by experts in adults as questionable results have been obtained in patients eating large meals very late. Fasting duration varies between 8 and 12 h in children and 3-4 h in infants less than 12 months old, depending on the clinical scenario.⁸¹

Drugs which influence gastric motor function should be avoided before the test. This includes established prokinetics as well as drugs with anticholinergic properties (e.g., tricyclic antidepressants), smooth muscle relaxants, and opioids. The duration of withdrawal depends on the half-life of the drug, 48–72 h are usually sufficient. However, in a patient with dyspeptic symptoms, who requires long-term medication with, for example, amitriptyline, it is not reasonable to alter the normal clinical situation by discontinuation of the drug before the test.

As ¹³C-GEBT are harmless, they can be performed repeatedly and have been used successfully to monitor drug effects in clinical studies and individual patients. ¹¹⁷⁻¹²⁰

How should test results be reported?

Recommendation 3.7 To allow for reliable interpretation of test results the following parameters should be reported: assessment of gastric emptying of solids versus liquids; marker substance; caloric content of the test meal; duration of breath sampling period; test result including normal values; and interpretation (accelerated, normal, delayed gastric emptying) with T½ being the best established and most widely used parameter (94%, 0%; Q:D S:B).

Results can be calculated from breath test data alone, not corrected for scintigraphy 15,20 or given as scintigraphic equivalent values (for T½) according to Ghoos et al. 13

Recommendation 3.8 The test report could be complemented by including composition and preparation of the test meal, breath sampling intervals, methods used for analysis, and additional

gastric emptying parameters such as T lag and GEC (gastric emptying coefficient) with normal values (86%, 0%; Q:D S:C).

These recommendations are based on the expert consensus that a minimum of information on methodology and test results including normal ranges and interpretation (accelerated/normal/delayed gastric emptying of solids/liquids) is required for reliable interpretation. Further information on methodological aspects and inclusion of several outcome parameters could improve diagnostic gain.

¹³C-PANCREATIC FUNCTION TESTS

In adults, chronic pancreatitis, pancreatic cancer, and surgical procedures are the most common causes of PEI, whereas in children, cystic fibrosis is of particular relevance. In these diseases, inflammatory destruction, pancreatic atrophy, ductal obstruction, or resection of pancreatic tissue lead to decreased exocrine secretion. Furthermore, PEI can be caused by regulatory imbalances in the presence of a normal pancreas, such as reduced hormonal and vagal stimulation of pancreatic secretion or inactivation of pancreatic enzymes in the intestinal lumen due to hyperchlorhydria (e.g., Zollinger-Ellison syndrome).⁴⁹

Due to the large reserve capacity of the pancreas, mild to moderate exocrine insufficiency is frequently not associated with clinical symptoms, and overt steatorrhea is not expected unless the secretion of pancreatic lipase is reduced to less than 10% of normal (severe PEI). However, patients with "compensated" PEI also have an increased risk of nutritional deficiencies, in particular, of lipid-soluble vitamins with respective clinical consequences.⁷

Indications for pancreatic function testing in adults and children as recommended by international guidelines are given in Table 7.7.121-124

When should ¹³C-PFBT be utilized?

Recommendation 4.1 The ¹³C-mixed triglyceride breath test is to be regarded as an established alternative to other non-invasive and invasive pancreatic function tests (97%, 0%; Q:A S:A).

Several breath tests using ¹³C-labeled lipids, proteins, or complex carbohydrates have been developed that indirectly assess pancreatic lipase, protease, and amylase activities. ^{11,12,14,100,127-136} Tests investigating lipid digestion and absorption are preferred because, from a clinical point of view, steatorrhea is by far the most important digestive malfunction in PEI: It is generally more severe and develops several years prior to malabsorption of protein or starch⁴⁹ and can be associated with decreased absorption of lipid-soluble vitamins, as mentioned above. Indeed, a ¹³C-breath test using naturally enriched maize starch lacks sensitivity and specificity for diagnosis of PEI. ¹³⁶

TABLE 7 Indications for pancreatic function testing in adults and children^{7,121-126}

Grade of recommendation	n	
Adults	High	Suspected pancreatic exocrine insufficiency in patients with pancreatic disease/after pancreatic resection
		Patients with chronic pancreatitis at the time of diagnosis (and annually thereafter if not tested positive for pancreatic exocrine insufficiency)
	Moderate	Monitoring of pancreatic enzyme replacement therapy in patients with an inadequate therapeutic response ^a
		Suspected pancreatic exocrine insufficiency without evidence of pancreatic disease
		Differential diagnosis of chronic diarrhea
Children	High	Screening for pancreatic exocrine insufficiency of children with chronic pancreatitis every 6–12 months
		Newly diagnosed CF
		Every 3–12 months (age dependent) in CF patients with pancreatic sufficiency at time of diagnosis
		In CF patient with subnormal weight development

Abbreviation: CF, cystic fibrosis.

Tests with various ¹³C-labeled lipids have been investigated, in which triolein, trioctanoin, tripalmitin, cholesteryl-octanoate, and mixed triglycerides (MTG) are used to generate marker substances. 128 Currently, the original 127 or modified versions 11,12,14,133 of the ¹³C-MTGBT developed by Vantrappen et al. are almost exclusively used in the clinic because of practical advantages (e.g., shorter labeled fatty acid allowing for shorter breath sampling period). It is based on the principle that intestinal triglyceride absorption requires prior hydrolysis by pancreatic lipase to produce free fatty acids and mono-acyl-glycerol. These metabolites are incorporated into micelles, absorbed, resynthesized, and transported to the liver. Hepatic enzymes subsequently release fatty acids, including ¹³C-octanoic acid, that is specifically bound to the Sn-2 position of ¹³C-mixed triglycerides (¹³C-MTG: 1,3 distearyl 2[¹³Coctanoyl] glycerol). ¹³C-octanoic acid then undergoes β-oxidation, which results in the formation of ¹³CO₂, which is absorbed into the bloodstream, transported to the lung, and exhaled. The increase in ¹³CO₂-concentration in breath thus correlates with pancreatic lipase secretion.

Direct comparison with the reference standard (determination of pancreatic enzyme and/or bicarbonate output in duodenal aspirates following exogenous stimulation with secretin \pm cerulein) demonstrates high sensitivity for severe exocrine insufficiency (90%–100%) with specificity ranging between 80% and 90% in adults. ^{121,137} A modified test using comparably high lipid loads in subjects explicitly avoiding physical activity during testing reached high sensitivity and specificity rates even in mild to moderate PEI (100% and 92%, respectively). ¹²

However, as evident from the test principle, the ¹³C-MTGBT is a test of lipid digestion and absorption. Therefore, the ¹³C-MTGBT is accepted as an appropriate alternative to the coefficient of fat

absorption, both for the diagnosis of PEI and for evaluating the efficacy of pancreatic enzyme replacement therapy (PERT) in clinical practice. 7.11,132,133,138 On the other hand, 13CO₂ exhalation is not only decreased by lipase deficiency but also by other causes of lipid malabsorption, for example, celiac disease, short bowel syndrome, or postcibal asynchrony following gastric resection. 139 Thus, the specificity of the test for the differential diagnosis of chronic diarrhea is limited.

Unfortunately, fecal elastase, the other noninvasive test which is predominantly used for pancreatic function testing in the clinic, also has limited specificity in these cases. 137 Direct comparisons between 13 C-MTGBT and fecal elastase favor the 13 C-MTGBT for diagnosis of steatorrhea, 140 and generally in patients with chronic pancreatitis after pancreatic resections. 141

In infants with cystic fibrosis, sensitivity of the ¹³C-MTGBT for diagnosis of steatorrhea was high, but specificity was low. ¹⁴² Thus, the test has been mainly, though not exclusively, ^{143,144} used to evaluate the efficacy of PERT in children with cystic fibrosis. ¹⁴⁵

How should ¹³C-PFBT be performed?

Recommendation 4.2 ¹³C-MTGBT performed for clinical reasons have to adhere strictly to standardized study protocols adequately validated in a representative patient population. This refers to preparation of the test meal as well as test performance and evaluation of test results (for examples compare Table 8) (100%, 0%; Q:C S:A).

Recommendation 4.3 Adult patients, adolescents and older children have to be fasted overnight (97%, 0%; Q:C S:A).

^aonly applicable for indirect tests measuring digestion/absorption.

TABLE 8 Established test protocols for ¹³C-MTGBT in adults and children

Remarks	29 pts with pancreatic disease, controls: 25 healthy subjects + 22 pts with nonpancreatic steatorrhea comparison with both, stimulated duodenal lipase output (reference standard for pancreatic secretion) and quantitative fecal fat (reference standard for steatorrhea) Effect of PERT demonstrated in subgroup of pancreatic pts	Developed using quantitative fecal fat (reference standard for steatorrhea) for comparison in healty volunteers (N = 10) and chronic pancreatitis patients with (N = 16) or without (N = 4) PEI, validated in 78 pts with advanced CP, also shown to be of value for monitoring of PERT efficacy, and to correlate with the nutritional status and the severity of chronic pancreatitis	Validated using secretin test (reference standard for pancreatic secretion) for comparison in HC and patients with pancreatic disease (N = 19), also detects mild and moderate PEI	Evaluated in 200 pts undergoing both, ¹³ C-MTGT and ¹³ C-GEBT. More convenient, but decreasing duration of the test associated with lower diagnostic accuracy. Tests with less than 4 h duration are markedly influenced by gastric emptying time
Validity	For detection of PEI (decreased lipase output): SENS 89%, SPEC 81% PPV 63% NPV 95% (control pts with nonpancreatic steatorrhea included)	SENS 93%, SPEC 92% ACC 92% (4h-test is associated with slightly lower diag- nostic ACC.: SENS 91%, SPEC 89%)	SENS 100%, SPEC 92% verses secretin-test	SENS 88% SPEC 94%, versus 6 h test version
Endpoints and normal values	Cumulative ¹³ C-recovery, normal (estimated from fig 4: Lowest value obtained in HC): >23% of dose	Cumulative ¹³ C-recovery, normal >29% of dose (>19% for 4 h test duration)	Cumulative ¹³ C-recovery, normal >26.8% of dose	Cumulative ¹³ C- recovery, normal >13.8% of dose
Breath sampling	At baseline and at 30 min intervals for 6 h pp	At baseline and at 30 min intervals for 6 h pp	At baseline and at 30 min intervals for 6 h pp	At baseline and at 30 min intervals for 4 h pp
Test meal	100 g of toast with 0.25 g of butter per kilogram of body weight, plus 16 mg ¹³ C -MTG per gram of butter	40 g bread, 20 g butter 200 ml water, ¹³ C-MTG spread on butter (total fat content 16 g) plus 250 mg ¹³ C-MTG 10 mg metoclopramide 20 min before meal ingestion	Two slices of white bread, 20 g butter, 30 g chocolate cream (31 g fat/100 g) carefully mixed with 250 mg ¹³ C-MTG (total fat content 26 g)	Two slices of white bread, 20 g butter, 30 g chocolate cream (31 g fat/100 g) carefully mixed with 250 mg ¹³ C-MTG (total fat content 26 g)
Age	Adults	Adults	Adults	Adults
	Vantrappen et al. 1989 ¹²⁷	Dominguez-Munoz et al. 2015 ¹³²	Keller et al. 2011 ¹²	Keller et al. 2014 ¹⁴

rable 8 (Continued)

/an

	Age	Test meal	Breath sampling	Endpoints and normal values	Validity	Remarks
n Dijk-van Aalst et al. 2001 ¹⁴⁶	12 premature infants, 12 full-term infants (1-6 months), 20 children (3-10 years), 20 teenagers (11-17 years)	12 premature infants, 12 Infants: Formula with low Two samples at baseline, Cumulative 6h- ¹³ CO ₂ -full-term infants (1- ¹³ C content (e.g., further samples at excretion (% of dos 6 months), 20 NAN1 (Nestlé), Pre-15 min intervals for administered) children (3- ¹⁰⁰ mg ¹³ C-MTG and true infants: 2.3.9 ± 5.2% administered) mean ± SD: Premateenagers (11- ¹⁰⁰ glycol 3350 > 3 years: slice of white bread with 5 g butter and 15 g chocolate paste, mixed with 250 mg ¹³ C-MTG, 100 ml whole-fat milk	Two samples at baseline, further samples at 15 min intervals for 6 h pp	Cumulative $6h^{-13}CO_2$ - excretion (% of dose administered) mean \pm SD: Premature infants: $23.9 \pm 5.2\%$ Full-term infants: $31.9 \pm 7.7\%$ Children: $32.5 \pm 5.3\%$ Teenagers: $28.0 \pm 5.4\%$		Mean value for healthy adults: 35.6%, lower limit of normal 22.8%

Abbreviations: ¹³C-GEBT, ¹³C-gastric emptying breath test; ¹³C-MTG, mixed triglycerides; ¹³C-MTGBT, ¹³C-mixed triglyceride breath test; ACC, accuracy; HC, healthy controls; PEI, pancreatic exocrine nsufficiency; PERT, pancreas enzyme replacement therapy; pp, postprandially; SENS, sensitivity; SPEC, specificity. **Recommendation 4.4** A shorter fasting period can be sufficient in very young children (95%, 0%; Q:C S:B).

Recommendation 4.5 Before and during the test, precautions as described in general methodology have to be taken into account (avoidance of ¹³C-rich food, avoidance of physical activity and ¹³C-rich infusions during test) (100%, 0%; Q:D S:C).

Compared to ¹³C-GEBT, ¹³C-MTGBT uses significantly lower amounts of ¹³C resulting in lower DOB values and lower cumulative ¹³C-exhalation rates. Consequently, avoidance of confounders such as ¹³C-rich food is of particular importance when conducting the ¹³C-MTGBT. Test results were not affected by liver disease in children with cystic fibrosis.¹⁴⁷

While adults and older children are required to fast overnight prior to breath testing, 11,12,127,132,133 a shorter period may be sufficient in very young children. In analogy to gastric emptying studies, 3-4 h of fasting are recommended in infants less than 12 months old, depending on the clinical scenario. 81

A mathematical correction for altered $\rm CO_2$ -production rates in physically active subjects using heart rate recordings improves diagnostic accuracy of the $^{13}\text{C-MTGBT}$ in non-resting subjects. 148 However, under clinical conditions, it is pragmatic to request that patients remain seated throughout the test.

Recommendation 4.6 Drugs with potential influence on gastrointestinal transit and/or lipid digestion and absorption should be avoided before the test. The duration of withdrawal depends on the half-life of the drug, 48–72 h are usually sufficient (91%, 3%; Q:C S:B).

Recommendation 4.7 Pancreatic enzyme replacement therapy should be avoided 2 days before the test, unless the test is performed for monitoring the therapeutic efficacy of pancreatic enzyme replacement therapy (97%, 0%; Q:B S:B).

Intraluminal lipolysis of the marker and lipid absorption require prior delivery of the test meal into the duodenum, adequate mixing of the marker, meal nutrients and digestive secretions and regulated transport of chyme through the small intestine. Therefore, drugs, which accelerate (e.g., prokinetics, potent laxatives) or markedly inhibit (e.g., opioids) GI transit, may alter test results and should be avoided before the test, unless essential long-term medication is concerned. Depending on the half-life of the drug, 48–72 h of withdrawal are usually sufficient. However, some test protocols deliberately include administration of metoclopramide in order to accelerate gastric emptying and to shorten the required breath sampling period. 11.132,149

If the ¹³C-MTGBT is performed to assess the presence of PEI, it is obvious that subjects have to stop oral PERT prior to testing. To account for potentially slow GI transit, it is recommended that patients abstain from PERT for 2 days before the test.

On the other hand, ¹³C-MTGBT is a suitable and established alternative to fecal fat measurement for monitoring the efficacy of

PERT. ^{11,132} To answer this question, subjects are allowed to continue with their usual enzyme replacement therapy. Indeed, it has been shown that PERT can be optimized by repetitive testing with increasing PERT doses until normal ¹³C-MTGBT results are achieved. By this, a significant increase of body weight was observed. ¹¹

How should test results be reported?

Recommendation 4.8 To allow for reliable interpretation of test results, the following parameters should be reported: marker substance, test results and interpretation including normal values. Cumulative ¹³C-recovery rate (in % of dose administered) represents the established main outcome parameter (94%, 0%; O:D S:B).

Recommendation 4.9 The test report can be complemented by including the ¹³C-exhalation curve, composition and preparation of the test meal, breath sampling intervals, methods used for analysis and interpretation of the ¹³C-exhalation with respect to the clinical context (91%, 0%; Q:D S:B).

These recommendations are based on the expert consensus that a minimum of information on test methodology and test results including normal ranges and interpretation (normal/decreased intestinal lipolysis compatible with PEI) is required for reliable interpretation. Further information regarding methodological aspects and interpretation of test results with respect to the individual clinical context could further improve the diagnostic validity.

¹³C-LIVER FUNCTION BREATH TESTS

Established parameters for the assessment of liver function under routine clinical conditions are measurement of bilirubin, albumin, liver enzymes and parameters of coagulation factor synthesis in serum or plasma, respectively. Clinical prognostic grading systems (e.g., Child-Pugh score, Model for End-stage Liver Disease score) combine several of these biochemical parameters including clinical symptoms of advanced liver cirrhosis. ¹⁵⁰

In contrast to these "static" liver function tests, "dynamic" quantitative tests measure the elimination of a substance, which is cleared and/or metabolized almost exclusively by the liver via specific metabolic pathways in subcellular compartments, for example, by cytochromes for microsomal liver function or cytosolic or mitochondrial enzymes. 151 Accordingly, these tests constitute a more accurate measure of the specific aspects of liver function. Established dynamic quantitative liver function tests are the indocyanine green clearance test and the galactose elimination capacity test. 150 13 C-LFBT also represent dynamic tests with oral consumption or intravenous application of the marker substance and measurement of the end product of hepatic metabolism, that is exhaled 13 CO₂. Some modifications have been shown to detect early changes in liver metabolic capacity in patients, prior to the

presence of structural damage to the liver (i.e., inflammation, fibrosis), 152 though with limited sensitivities and specificities. $^{153-160}$ Still, the dynamic nature of 13 C-LFBT, their possible versatility in terms of assessing a range of different liver functions, and the ease with which they can be repeated to follow relative changes in liver function with time, generally imply a marked potential for clinical application. 40

When should ¹³C-LFBT be utilized?

Recommendation 5.1 ¹³C-liver function breath tests could be used for measurement of various aspects of liver function in adults (94%, 0%; Q:A S:C).

Recommendation 5.2 Presently, due to very limited evidence, performance of ¹³C-liver function breath tests for clinical reasons cannot be recommended in children (95%. 0%: O:D S:B)

Different ¹³C-LFBT have been developed for assessment of hepatic mitochondrial (substrates: ¹³C-ketoisocaproate, [methyl-¹³C]methionine), microsomal (13C-methacetin, 13C-aminopyrine, [3-methyl-13C]-caffeine) and cytosolic (13C-phenylalanine) function. 152,157,159-¹⁶⁸ These have been mainly used in patients with liver fibrosis and cirrhosis due to nonalcoholic fatty liver disease, 152,157,159-^{161,163,164,169,170} chronic hepatitis C, ^{162,165} cirrhosis, ^{166,171} and hepatocellular carcinoma. 172 Additional potential fields of application are in steatohepatitis. 156,157,159 fatty liver. 160,163,169 and assessment of prognosis in chronic liver diseases, in general: In a 7-year prospective follow-up study in 132 patients with chronic HCV infection, the 4'-O-13C-methacetin breath was not inferior to liver biopsy in predicting liver-related death and transplantation. 173 Although several studies have reported close correlations between ¹³C-LFBT and histological alterations or other established parameters in adults, they have still not entered the mainstream of clinical practice but are used exclusively by highly specialized centers. For instance, the effect of transarterial chemoembolization on liver function has been monitored by ¹³C-methacetin test, ¹⁷⁴ and several studies suggest that this test could also be used for planning of hepatic resections. 175-177

Given that, as yet, few studies have been performed involving small groups of children with rare diseases, $^{178-181}$ the use of 13 C-LFBT for clinical purposes cannot be recommended in the pediatric population, so far.

How should ¹³C-LFBT be performed?

Recommendation 5.3 ¹³C-liver function breath tests performed for clinical purposes have to adhere strictly to standardized and adequately validated study protocols. This refers to preparation of the patient, test meal/solution, test performance and evaluation of test results (for examples compare Table 9) (100%, 0%; Q:C S:A).

TABLE 9 Examples of validated test protocols for ¹³C-liver function breath tests in adults

	Estimated parameter	Marker and test solution	Breath sampling	Endpoints and normal values	Validity	Remarks
Afolabi et al 2018 ¹⁶¹	Hepatic mitochondrial function	1 mg/kg body weight of ¹³ C-ketoisocaproate plus 20 mg/kg body weight L-leucine dissolved in 200 ml of water	At baseline and at 10 min intervals for 60 min pp	Cumulative ¹³ C- recovery, normal >21% of dose		Validated in 11 HC and 77 pts with NAFLD, SENS and SPEC to detect significant fibrosis was not determined
Portincasa et al 2006 ¹⁵⁷	Hepatic mitochondrial function	1 mg/kg body weight of 13C-ketoisocaproate plus 1 g of L-leucine dissolved in 200 ml of water	At baseline and at 10 min intervals postprandially for 60 min pp	Cumulative ¹³ C- recovery, normal >14% of dose	Diagnostic accuracy at identifying pts with NASH (cut-off value 9.6%):SENS 68%, SPEC 94%,PPV 90%, NPV 73%	Validated in 28 HC and 39 pts with NAFLDThe test was also able to discriminate fibrosis stages in patients with NASH
Banasch et al 2011 ¹⁵⁶	Hepatic mitochondrial function	2 mg/kg body weight [methyl- ¹³ C]- methionine dissolved in 100 ml of water, prior oral consumption of 200 ml of orange juice	At baseline and at 10 min intervals for 90 min pp	Cumulative ¹³ C-recovery, normal >6.1% of dose	Cut-off value <4.20% of dose for separation of pts with NASH from non-NASH, SENS 81%; SPEC 76%	Validated in 118 pts with NAFLD and 18 HCTest predicts higher stages of disease activity
Korkmaz et al. 2015 ¹⁶⁰	Hepatic mitochondrial function	2 mg/kg body weight [methyl- ¹³ C]- methionine dissolved in 100 ml water	At baseline and at 10 min intervals for 90 min pp	Cumulative ¹³ C-recovery, normal >6.2% of dose	Cut off value <3.71%: SENS 95% SPEC 88% for differentiating advanced liver fibrosis (F2-3) from mild (F0-1) fibrosis	Validated in 164 pts with NAFLD and 56 HC
Fierbinteanu-Braticevici et al.2013 ¹⁶³	Hepatic microsomal function	Fixed dose of 75 mg ¹³ C- methacetin dissolved in 200 ml of water	At baseline and at 10 min intervals for 60 min pp	Cumulative ¹³ C- recovery, normal >22% of dose	Cut off value <15.2%, SENS 91%, SPEC 82% at detecting significant fibrosis $(F \ge 2)$	Validated in 90 pts with NAFLD and 20 HC
Park et al 2003 ¹⁶⁷	Hepatic microsomal function	2 mg/kg body weight of [3-methyl- ¹³ C]-caffeine dissolved in 30 ml of water, followed by 40 ml water wash of the container	At baseline and at 10 min intervals for 60 min pp	Cumulative ¹³ C- recovery, normal >2.3% of dose	Cut off value <1.85%: 79% SENS 80% SPEC for detecting NASH	Validated in 48 pts with NAFLD, 48 patients with chronic hepatitis B and 24 HC subjects. Results reflect the extent of hepatic fibrosis

Abbreviations: HC, healthy controls; NAFLD, non-alcoholic Fatty Liver disease; NASH, nonalcoholic steatohepatitis; NPV, negative predictive value; PPV, positive predictive value; pp, postprandial; SENS, sensitivity, SPEC, specificity.

For most validated test protocols, the marker substance is dissolved in up to 200 ml of water or unsweetened tea; breath samples are collected before drinking of the test solution and at regular intervals for up to 1.5 h postprandially. As shown in Table 9, cumulative percentage of the marker dose exhaled at the end of the observation period is the primary outcome parameter of the vast majority of tests. The list of tests detailed in the table is incomplete but concentrates on examples for test protocols using different marker substances, which have been validated in an adequate patient population and have established normal or cut off values for specific diagnoses. Apart from tests with oral marker application, the LIMAx test¹⁷⁷ is commercially available in some European countries, which uses intravenous application of 4'-O-¹³C-methacetin. For this test, the manufacturer's instructions have to be followed exactly.

Recommendation 5.4 Before and during the test precautions as described in general methodology have to be observed (avoidance of ¹³C-rich food, of physical activity, and of ¹³C-rich infusions) (100%, 0%; Q:C S:A).

Recommendation 5.5 Drugs with potential influence on gastrointestinal transit and/or cytochrome p450 metabolism should be avoided before and during the test (96%, 0%; Q:C S:B).

Recommendation 5.6 As far as essential long-term medication is concerned, it can be necessary to perform ¹³C-liver function breath tests despite ongoing medication as long as these medications are not administered during the test period (93%, 0%; Q:B S:B).

As discussed before, potential confounders of test results detailed in general methodology have to be avoided. For $^{13}\text{C-LFBT}$ that investigate cytochrome P₄₅₀ enzymes, drugs with potential influence on cytochrome P₄₅₀ metabolism are of particular relevance, $^{182-184}$ and it may be necessary to perform a test despite ongoing medication.

Especially for the most frequently applied test substrate ¹³C-methacetin, strong CYP1A2 inhibitors (e.g., ciprofloxacin, fluvoxamin) will have an influence on hepatic methacetin metabolism, although no interaction studies have been systematically performed. Influences of weak CYP1A2 inhibitors (e.g., norfloxacin, propranolol) have not been reported. However, due to the rather short duration of the tests (no longer than 2 h for most protocols), administration during the test period can usually be circumvented, even if essential long-term medication is concerned.

How should test results be reported?

Recommendation 5.7 To allow for reliable interpretation of test results the following parameters should be reported: marker substance; test results and interpretation including normal values (93%, 0%; Q:D S:B).

Recommendation 5.8 The test report could be complemented by including the ¹³C-exhalation curve, composition and preparation of the test meal/solution, breath sampling intervals, methods used for analysis, and interpretation of test results with respect to the clinical context (82%, 0%; Q:D S:C).

As for other ¹³C-breath tests, these recommendations are based on the expert consensus that a minimum of information on methodology and test results including normal ranges and interpretation is required for reliable interpretation. Further information regarding methodological aspects and interpretation of test results with respect to the individual clinical context could further improve the diagnostic validity.

OTHER ¹³C-BREATH TESTS

The ¹³C-breath tests discussed in the previous chapters have gastroenterological indications and are either clinically established or at least used regularly by specialized centers. Apart from these, there is a multitude of other test options at various developmental stages for gastroenterological indications and ¹³C-breath tests, which are used by other medical specialties. Examples for these are detailed in Table 10 and demonstrate the very broad potential applicability of ¹³C-breath test technology.

CONCLUSIONS AND FUTURE PERSPECTIVES

This consensus-based clinical practice guideline aims to assist physicians and provide them with the information required to perform high quality ¹³C-breath tests for patients with various GI symptoms and diseases. In a consensus process, in which representatives from all regions of Europe and specialists representing European scientific societies participated, the available evidence was evaluated, taking account of local facilities, diverse clinical practice, and health care environments. Patient involvement was not covered in our guideline because of its focus on diagnostic recommendations. However, after publication of the guideline, a qualitative and quantitative assessment of guideline adoption into clinical practice is planned, including patient and public involvement.

The guideline gives an overview over general methodology of 13 C-breath testing and provides recommendations for the use of 13 C-breath tests to diagnose *H. pylori* infection, measure gastric emptying time and monitor pancreatic exocrine and liver function in adult and pediatric patients. Other potential applications of 13 C-breath testing are summarized briefly. The recommendations specifically detail when and how individual 13 C-breath tests should be performed including examples for well-established test protocols.

¹³C-breath tests are indirect tests and mostly require several digestive and metabolic steps that ultimately lead to the exhalation of ¹³C-enriched breath, including intestinal absorption, hepatic metabolization, and pulmonary excretion. Therefore, demographic

TABLE 10 Examples of other ¹³C-breath tests

References	Estimated parameter	Remarks
(185-194)	Insulin resistance, metabolic syndrome	$^{13}\text{C-glucose}$ + unmarked glucose (e.g., 75 g) dissolved in water, outcome parameters: Cumulative $^{13}\text{CO}_2\text{-exhalation}$, maximal DOB-value, $^{13}\text{C-exhalation}$ at specific time point
		¹³ C-exhalation is significantly decreased in patients with impaired insulin resistance/ metabolic syndrome, site of labeled carbon atom in glucose molecule is relevant for metabolic results
		Breath test described as valid surrogate index of clamp-derived measures of insulin resistance, with good accuracy and precision, but controversial findings regarding applicability as screening tool (insufficient accuracy) in minority of studies
		Non-invasive tool for early diagnosis and follow up of patients in high-risk groups
(195-198)	SIBO, OCTT	Lactose-[¹³ C]ureide, moderate sensitivity (66.7%), high specificity (100%) for SIBO in comparison with culture from jejunal aspirates, glucose-[¹³ C]ureide could be used alternatively
		Bacterial metabolism of marker is prerequisite for $^{13}\mathrm{CO}_2$ -exhalation, increase caused by SIBO or arrival of the marker substance in the colon
		Predosing with unlabelled ureide is frequently part of the protocol, induces enzyme induction, increases rate of $^{13}\mathrm{CO}_2$ -recovery but complicates procedure
(136,199-201)	PEI (amylase deficiency) Congenital sucrase-isomaltase deficiency	Maize starch is naturally enriched with 13 C, 50–100 g used as pancreatic function test, but limited sensitivity (up to 77%) and specificity (up to 74%) in comparison with "gold standard" (secretin cerulein test)
		$^{13}\mbox{CO}_2$ release following labelled starch or sucrose is decreased in patients with sucrase-isomaltase deficiency
(202-210)	Activity of drug metabolizing enzymes	13 C-erythromycin applied iv (or orally) \rightarrow hepatic (and intestinal) cytochrome P450 3A4 and 3A5 (CYP3A4/5) activity in vivo, pilot studies
		$^{13}\text{C-pantoprazole} \rightarrow \text{CYP2C19}$ activity, can predict the anti-platelet efficacy of clopidogrel (which is metabolized by CYP2C19 to its active metabolite) and high metabolization rate of PPI
		¹³ C-dextromethorphan → CYP2D6 activity (e.g., relevant for activation of the prodrug tamoxifen), [2-(13)C]uracil → dihydropyrimidine dehydrogenase (DPD) deficiency to predict 5-fluorouracil dose-related toxicity

Abbreviations: DOB, δ over baseline; OCTT, oro-cecal transit time; PEI, pancreatic exocrine insufficiency; PPI, proton pump inhibitor; SIBO, small intestinal bacterial overgrowth.

parameters, concomitant diseases, and/or medication affecting these processes may confound breath test results. In fact, studies have shown that the influence of concomitant diseases on many ¹³C-breath tests is small or can be prevented by adaptation of the breath test protocol. Information on general ¹³C-breath test methodology as detailed in the guideline will help clinicians to optimize patient selection and preparation and to better understand and interpret breath test results for individual patients.

As depicted in Figure 2, ¹³C-UBT has excellent sensitivity and specificity for diagnosis of *H. pylori* infection, superior to most other tests, it is highly standardized and well established in the clinic. ¹³C-GEBT and ¹³C-PFBT are validated and accepted tests for diagnosis of gastric emptying disturbances and PEI, respectively. So far, their wide spread use has been hampered by the lack of standardization among centers, which we aim to improve by this guideline. Moreover, the tests are currently not commercially available in Europe, so that test meals need to be prepared on site. ¹³C-LFBT are currently used by highly specialized centers, only. They have the potential to measure various

aspects of liver function, including assessment of prognosis in chronic liver disease. However, their clinical role, in general, and applicability in pediatric patients, in particular, still need to be established.

Apart from identification of this and other areas of future research, the guideline should improve pan-European harmonization of diagnostic approaches to symptoms and disorders, which are very common in specialist and primary care gastroenterology practice, both in adult and pediatric patients.

ACKNOWLEDGMENTS

This guideline was developed with the support of an UEG Activity Grant.

CONFLICT OF INTEREST

Oliver Goetze received financial support/honoraria for clinical studies/lectures from Kibion, Mayoly Spindler Laboratories. Stephan L. Haas received honoraria by Mylan for oral presentations. The other authors have nothing to disclose.

13CO₂ ¹³C-LFBT: → Liver function ¹³C-UBT: Helicobacter pylori infection Substrates for estimation of hepatic mitochondrial 75 mg ¹³C-urea, dissolved in citric acid solution (e.g. ¹³C-ketoisocaproate) or microsomal · Test kits commercially available (e.g. ¹³C-methacetin) function · 2 point measurement: baseline and at 15 min · Most tests: oral application of substrate in (children) or 30 min (adults) water/tea; i.v. ¹³C-methacetin test commercially • ≥95% sensitivity and specificity in adults and children available in some European countries Highest accuracy among non-invasive H. pylori tests, • Prolonged measurements over 60-90 min more accurate than rapid urease test with collection of samples at ~10 min intervals Results correlate with histology in adults, but tests only used by highly specialized centers. Not recommended in children because of lack of data. ¹³C-PFBT: → Pancreatic exocrine function ³C-GEBT: Gastric emptying velocity 13CO₂ Egg-based test meals labeled with 13C-octanoic acid 13C-mixed triglycerides mixed with butter/chocolate or ¹³C-Spirulina platensis for ⊠gastric emptying of cream and spread on bread (formula in infants) · No standardized test kit commercially available solids, 13C-acetate preferred for gastric emptying of ³C-substrate adherence to validated study protocols liauids • In Europe, no standardized test kit commercially including established analysis methods available adherence to validated study protocols • Prolonged measurements over ~4-6 h with 13CO2 including established analysis methods required! collection of samples at (15-) 30 min intervals Prolonged measurements over ~2 h (liquids) to ~4 h Accepted pancreatic function test in adults (solids) with collection of samples at ~15 min intervals for solids, 5-10 min for early liquid gastric emptying and children, can also be used to monitor 13C-urea enzyme replacement therapy Accepted alternative to gastric emptyin scintigraphy in adults and children

FIGURE 2 General principle, performance and clinical role of 13 C-breath tests used in gastroenterology. 13 C-urea used for detection of 13 C-urea used for detection of 13 C-UBT, marked in green) is metabolized by bacterial urease to produce 13 CO₂, which is absorbed, transported to the lung (broken green arrow) and exhaled. For 13 C-gastric emptying breath tests (13 C-GEBT, marked in red), 13 C-pancreatic function breath test (13 C-PFBT, marked in blue), orally applied substrates or their metabolites are absorbed in the small intestine. Subsequently, they are transported to the liver where they undergo further metabolization with production of 13 CO₂, which is transported to the lung and exhaled

DATA AVAILABILITY STATEMENT

Data sharing is not applicable to this article as no new data were created or analyzed in this study.

ORCID

Jutta Keller https://orcid.org/0000-0002-5884-1115

Heinz F. Hammer https://orcid.org/0000-0001-8683-4951

Mark R. Fox https://orcid.org/0000-0003-4394-5584

REFERENCES

- Sonyi M, Keller J, Fox M, Hammer HF. Development of a multinational clinical practice guideline: a practical structured procedure. Basel: Digestive diseases; 2020.
- Keller J, Franke A, Storr M, Wiedbrauck F, Schirra J. Clinically relevant breath tests in gastroenterological diagnostics--recommendations of the German Society for Neurogastroenterology and Motility as well as the German Society for Digestive and Metabolic Diseases. Z Gastroenterol. 2005;43 (9):1071-90.
- Schoeller DA, Schneider JF, Solomons NW, Watkins JB, Klein PD. Clinical diagnosis with the stable isotope ¹³C in CO2 breath tests: methodology and fundamental considerations. J Lab Clin Med. 1977;90 (3):412-21.
- Malfertheiner P, Megraud F, O'Morain CA, Gisbert JP, Kuipers EJ, Axon AT, et al. Management of Helicobacter pylori infection-the Maastricht V/Florence consensus report. Gut. 2017;66 (1):6–30.

- Jones NL, Koletzko S, Goodman K, Bontems P, Cadranel S, Casswall T, et al. Joint ESPGHAN/NASPGHAN guidelines for the management of Helicobacter pylori in children and adolescents (update 2016). J Pediatr Gastroenterol Nutr. 2017;64 (6):991–1003.
- Keller J, Bassotti G, Clarke J, Dinning P, Fox M, Grover M, et al. Expert consensus document: advances in the diagnosis and classification of gastric and intestinal motility disorders. Nat Rev Gastroenterol Hepatol. 2018;15 (5):291–308.
- Lohr JM, Dominguez-Munoz E, Rosendahl J, Besselink M, Mayerle J, Lerch MM, et al. United European Gastroenterology evidence-based guidelines for the diagnosis and therapy of chronic pancreatitis (HaPanEU). United European Gastroenterol J. 2017;5 (2):153–99.
- Jones PJ, Leatherdale ST. Stable isotopes in clinical research: safety reaffirmed. Clin Sci (Lond). 1991;80 (4):277–80.
- Klein PD, Klein ER. Stable isotopes: origins and safety. J Clin Pharmacol. 1986;26 (6):378–82.
- Gregg CT, Hutson JY, Prine JR, Ott DG, Furchner JE. Substantial replacement of mammalian body carbon with carbon-13. Life Sci. 1973;13 (7):775–82.
- Dominguez-Munoz JE, Iglesias-Garcia J, Vilarino-Insua M, Iglesias-Rey M. ¹³C-mixed triglyceride breath test to assess oral enzyme substitution therapy in patients with chronic pancreatitis. Clin Gastroenterol Hepatol. 2007;5 (4):484-8.
- Keller J, Bruckel S, Jahr C, Layer P. A modified (1)(3)C-mixed triglyceride breath test detects moderate pancreatic exocrine insufficiency. Pancreas. 2011;40 (8):1201-5.
- Jonderko K, Dus Z, Szymszal M, Kasicka-Jonderko A, Blonska-Fajfrowska B. Normative values for the ¹³C-mixed triglyceride

breath test in two age groups. Med Sci Mon Int Med J Exp Clin Res. 2009:15 (5): Cr255-9.

- Keller J, Meier V, Wolfram KU, Rosien U, Layer P. Sensitivity and specificity of an abbreviated (13)C-mixed triglyceride breath test for measurement of pancreatic exocrine function. United European Gastroenterol J. 2014;2 (4):288-94.
- Leodolter A, Domínguez-Muñoz JE, von Arnim U, Kahl S, Peitz U, Malfertheiner P. Validity of a modified ¹³C-urea breath test for pre- and posttreatment diagnosis of Helicobacter pylori infection in the routine clinical setting. Am J Gastroenterol. 1999;94 (8):2100-4.
- Savarino V, Vigneri S, Celle G. The ¹³C urea breath test in the diagnosis of Helicobacter pylori infection. Gut. 1999;45 (Suppl 1): 118–22.
- Savarino V, Mela GS, Zentilin P, Bisso G, Pivari M, Mansi C, et al. Comparison of isotope ratio mass spectrometry and nondispersive isotope-selective infrared spectroscopy for 13C-urea breath test. Am J Gastroenterol. 1999:94 (5):1203-8.
- Mana F, Franken PR, Ham HR, Reynaert H, Urbain D. ¹³C urea breath test with nondispersive isotope-selective infrared spectrometry: reproducibility and importance of the fasting status. Helicobacter. 2000;5 (2):104–8.
- Gisbert JP, Gomollon F, Dominguez-Munoz JE, Borda F, Jimenez I, Vazquez MA, et al. [Comparison between two ¹³C-urea breath tests for the diagnosis of Helicobacter pylori infection: isotope ratio mass spectrometer versus infrared spectrometer]. Gastroenterol Hepatol. 2003;26 (3):141-6.
- Plavnik RG, Nevmerzhitsky VI, Butorova LI, Plavnik TE. [C assessment of mass-spectrometry and infrared spectrometry used in ¹³C-urea breath test for helicobacter pylori]. Klin Med (Mosk). 2015;93 (9):42–5.
- Kato M, Saito M, Fukuda S, Kato C, Ohara S, Hamada S, et al. ¹³C-Urea breath test, using a new compact nondispersive isotope-selective infrared spectrophotometer: comparison with mass spectrometry. J Gastroenterol. 2004;39 (7):629–34.
- Ferwana M, Abdulmajeed I, Alhajiahmed A, Madani W, Firwana B, Hasan R, et al. Accuracy of urea breath test in Helicobacter pylori infection: meta-analysis. World J Gastroenterol. 2015;21 (4):1305-14.
- Schadewaldt P, Schommartz B, Wienrich G, Brosicke H, Piolot R, Ziegler D. Application of isotope-selective nondispersive infrared spectrometry (IRIS) for evaluation of [¹³C]octanoic acid gastricemptying breath tests: comparison with isotope ratio-mass spectrometry (IRMS). Clin Chem. 1997;43 (3):518–22.
- 24. Goetze O, Selzner N, Fruehauf H, Fried M, Gerlach T, Mullhaupt B. ¹³C-methacetin breath test as a quantitative liver function test in patients with chronic hepatitis C infection: continuous automatic molecular correlation spectroscopy compared to isotopic ratio mass spectrometry. Aliment Pharmacol Ther. 2007;26 (2):305–11.
- 25. Boedeker C, Goetze O, Pfaffenbach B, Luypaerts A, Geypens B, Adamek RJ. ¹³C mixed-triglyceride breath test: isotope selective non-dispersive infrared spectrometry in comparison with isotope ratio mass spectrometry in volunteers and patients with chronic pancreatitis. Scand J Gastroenterol. 1999;34 (11):1153–6.
- Craig H. Isotopic standards for carbon and oxygen and correction factors for mass spectrometric analysis of carbon dioxide Geochemica. Cosmochimica Acta. 1957;12:181–6.
- van derMerwe N. Carbon isotopes, photosynthesis, and archaeology. Am Sci. 1982;70:596-606.
- Colaiocco Ferrante L, Papponetti M, Marcuccitti J, Neri M, Festi D.
 ¹³C-urea breath test for helicobacter pylori infection: stability of samples over time. Scand J Gastroenterol. 1999;34 (9):942–3.
- Irving CS, Wong WW, Shulman RJ, Smith EO, Klein PD. [¹³C]bicarbonate kinetics in humans: intra- vs. interindividual variations. Am J Physiol. 1983;245 (2):R190–202.

- Kindermann A, Demmelmair H, Koletzko B, Krauss-Etschmann S, Wiebecke B, Koletzko S. Influence of age on ¹³C-urea breath test results in children. J Pediatr Gastroenterol Nutr. 2000;30 (1):85–91.
- Caubet MS, Laplante A, Caille J, Brazier JL. [¹³C]aminopyrine and [¹³C]caffeine breath test: influence of gender, cigarette smoking and oral contraceptives intake. Isot Environ Health Stud. 2002;38 (2):71–7.
- Hellmig S, Von Schoning F, Gadow C, Katsoulis S, Hedderich J, Folsch UR, et al. Gastric emptying time of fluids and solids in healthy subjects determined by ¹³C breath tests: influence of age, sex and body mass index. J Gastroenterol Hepatol. 2006;21 (12):1832-8.
- Goetze O, Fox M, Kwiatek MA, Treier R, Schwizer W, Thumshirn M, et al. Effects of postgastric ¹³C-acetate processing on measurement of gastric emptying: a systematic investigation in health. Neuro Gastroenterol Motil. 2009;21 (10):1047-e85.
- Sonko BJ, Prentice AM, Coward WA, Murgatroyd PR, Goldberg GR.
 Dose-response relationship between fat ingestion and oxidation: quantitative estimation using whole-body calorimetry and ¹³C isotope ratio mass spectrometry. Eur J Clin Nutr. 2001;55 (1):10–8.
- McCue MD, Passement CA, Rodriguez M. The magnitude of the naturally occurring isotopic enrichment of ¹³C in exhaled CO2 is directly proportional to exercise intensity in humans. Comp Biochem Physiol Mol Integr Physiol. 2015;179:164–71.
- Keller J, Fliegner-Baia M, Layer P. Physical activity alters normal values of the "European standard" ¹³C-octanoic acid breath test. Gut. 2002;51 (Suppl III):A136.
- Bammens B, Evenepoel P, Rutgeerts P, Ghoos Y, Maes B. ¹³C-breath tests in peritoneal dialysis patients: influence of dialysis fluids. Eur J Gastroenterol Hepatol. 2003;15 (8):931–2.
- Eaton S, Pacilli M, Wood J, McHoney M, Corizia L, Kingsley C, et al. Factors affecting ¹³C-natural abundance measurement of breath carbon dioxide during surgery: absorption of carbon dioxide during endoscopic procedures. Rapid Commun Mass Spectrom. 2008;22 (11):1759–62.
- Keller J, Andresen V, Wolter J, Layer P, Camilleri M. Influence of clinical parameters on the results of ¹³C-octanoic acid breath tests: examination of different mathematical models in a large patient cohort. Neuro Gastroenterol Motil. 2009;21 (10):1039e83.
- Afolabi P, Wright M, Wootton SA, Jackson AA. Clinical utility of ¹³C-liver-function breath tests for assessment of hepatic function. Dig Dis Sci. 2013;58 (1):33–41.
- Westerterp KR, Plasqui G. Physical activity and human energy expenditure. Curr Opin Clin Nutr Metab Care. 2004;7 (6):607–13.
- Keim NL, Blanton CA, Kretsch MJ. America's obesity epidemic: measuring physical activity to promote an active lifestyle. J Am Diet Assoc. 2004;104 (9):1398–409.
- Costa RJS, Snipe RMJ, Kitic CM, Gibson PR. Systematic review: exercise-induced gastrointestinal syndrome-implications for health and intestinal disease. Aliment Pharmacol Ther. 2017;46 (3):246–65.
- van de Casteele M, Luypaerts A, Geypens B, Fevery J, Ghoos Y, Nevens F. Oxidative breakdown of octanoic acid is maintained in patients with cirrhosis despite advanced disease. Neuro Gastroenterol Motil. 2003;15 (2):113–20.
- 45. Mawatari H, Inamori M, Fujita K, Yoneda M, Iida H, Endo H, et al. The continuous real-time ¹³C-octanoate breath test for patients with nonalcoholic steatohepatitis using the BreathID system. Hepato-gastroenterology. 2009;56 (94-95):1436-8.
- Schneider AR, Kraut C, Lindenthal B, Braden B, Caspary WF, Stein J. Total body metabolism of ¹³C-octanoic acid is preserved in patients with non-alcoholic steatohepatitis, but differs between women and men. Eur J Gastroenterol Hepatol. 2005;17 (11):1181–4.

- Chapman MJ, Besanko LK, Burgstad CM, Fraser RJ, Bellon M, O'Connor S, et al. Gastric emptying of a liquid nutrient meal in the critically ill: relationship between scintigraphic and carbon breath test measurement. Gut. 2011;60 (10):1336-43.
- Maes BD, Ghoos YF, Geypens BJ, Hiele MI, Rutgeerts PJ. Relation between gastric emptying rate and rate of intraluminal lipolysis. Gut. 1996;38 (1):23-7.
- 49. Keller J, Layer P. Human pancreatic exocrine response to nutrients in health and disease. Gut. 2005;54 (Suppl 6):1–28.
- Best LM, Takwoingi Y, Siddique S, Selladurai A, Gandhi A, Low B, et al. Non-invasive diagnostic tests for Helicobacter pylori infection. Cochrane Database Syst Rev. 2018;3 (3):Cd012080.
- Chey WD, Wong BC. American College of Gastroenterology guideline on the management of Helicobacter pylori infection. Am J Gastroenterol. 2007;102 (8):1808–25.
- Leal YA, Flores LL, Fuentes-Panana EM, Cedillo-Rivera R, Torres J.
 ¹³C-urea breath test for the diagnosis of Helicobacter pylori infection in children: a systematic review and meta-analysis. Helicobacter. 2011;16 (4):327–37.
- Fischbach W, Malfertheiner P, Lynen Jansen P, Bolten W, Bornschein J, Buderus S, et al. S2k-Leitlinie Helicobacter pylori und gastroduodenale Ulkuskrankheit. Z Gastroenterol. 2017;54 (4):327–63.
- Leodolter A, Domínguez-Muñoz JE, Von Arnim U, Malfertheiner P. Citric acid or orange juice for the ¹³C-urea breath test: the impact of pH and gastric emptying. Aliment Pharmacol Ther. 1999;13 (8):1057-62.
- Domínguez-Muñoz JE, Leodolter A, Sauerbruch T, Malfertheiner P. A citric acid solution is an optimal test drink in the ¹³C-urea breath test for the diagnosis of Helicobacter pylori infection. Gut. 1997;40 (4):459-62.
- Elitsur Y, Tolia V, Gilger MA, Reeves-Garcia J, Schmidt-Sommerfeld E, Opekun AR, et al. Urea breath test in children: the United States prospective, multicenter study. Helicobacter. 2009:14 (2):134–40.
- Graham DY, Klein PD, Evans DJ Jr., Evans DG, Alpert LC, Opekun AR, et al. Campylobacter pylori detected noninvasively by the ¹³Curea breath test. Lancet (London, England). 1987;1 (8543):1174-7.
- Di Rienzo TA, D'Angelo G, Ojetti V, Campanale MC, Tortora A, Cesario V, et al. ¹³C-Urea breath test for the diagnosis of Helicobacter pylori infection. Eur Rev Med Pharmacol Sci. 2013;17 (Suppl 2):51–8.
- Logan RP, Polson RJ, Misiewicz JJ, Rao G, Karim NQ, Newell D, et al. Simplified single sample ¹³Carbon urea breath test for Helicobacter pylori: comparison with histology, culture, and ELISA serology. Gut. 1991;32 (12):1461–4.
- Gisbert JP, Pajares JM. Review article: ¹³C-urea breath test in the diagnosis of Helicobacter pylori infection -- a critical review. Aliment Pharmacol Ther. 2004;20 (10):1001–17.
- Rowland M, Lambert I, Gormally S, Daly LE, Thomas JE, Hetherington C, et al. Carbon 13-labeled urea breath test for the diagnosis of Helicobacter pylori infection in children. J Pediatr. 1997;131 (6):815–20.
- Gatta L, Vakil N, Ricci C, Osborn JF, Tampieri A, Perna F, et al. Effect of proton pump inhibitors and antacid therapy on ¹³C urea breath tests and stool test for Helicobacter pylori infection. Am J Gastroenterol. 2004;99 (5):823–9.
- 63. Graham DY, Opekun AR, Hammoud F, Yamaoka Y, Reddy R, Osato MS, et al. Studies regarding the mechanism of false negative urea breath tests with proton pump inhibitors. Am J Gastroenterol. 2003;98 (5):1005–9.
- Savarino V, Tracci D, Dulbecco P, Mele MR, Zentilin P, Mansi C, et al. Negative effect of ranitidine on the results of urea breath test for the diagnosis of Helicobacter pylori. Am J Gastroenterol. 2001;96 (2):348–52.

- Camilleri M, Parkman HP, Shafi MA, Abell TL, Gerson L. Clinical guideline: management of gastroparesis. Am J Gastroenterol. 2013;108 (1):18–37.
- Parkman HP, Camilleri M, Farrugia G, McCallum RW, Bharucha AE, Mayer EA, et al. Gastroparesis and functional dyspepsia: excerpts from the AGA/ANMS meeting. Neuro Gastroenterol Motil. 2010;22 (2):113–33.
- Vijayvargiya P, Jameie-Oskooei S, Camilleri M, Chedid V, Erwin PJ, Murad MH. Association between delayed gastric emptying and upper gastrointestinal symptoms: a systematic review and metaanalysis. Gut. 2019;68 (5):804–13.
- Rodriguez L, Irani K, Jiang H, Goldstein AM. Clinical presentation, response to therapy, and outcome of gastroparesis in children. J Pediatr Gastroenterol Nutr. 2012;55 (2):185–90.
- Waseem S, Islam S, Kahn G, Moshiree B, Talley NJ. Spectrum of gastroparesis in children. J Pediatr Gastroenterol Nutr. 2012;55 (2):166-72.
- Saliakellis E, Fotoulaki M. Gastroparesis in children. Ann Gastroenterol. 2013;26 (3):204-11.
- Kempler P, Amarenco G, Freeman R, Frontoni S, Horowitz M, Stevens M, et al. Management strategies for gastrointestinal, erectile, bladder, and sudomotor dysfunction in patients with diabetes. Diabetes Metab Res Rev. 2011;27 (7):665–77.
- Kovacic K, Elfar W, Rosen JM, Yacob D, Raynor J, Mostamand S, et al. Update on pediatric gastroparesis: a review of the published literature and recommendations for future research. Neuro Gastroenterol Motil. 2020;32 (3): e13780.
- Abell TL, Camilleri M, Donohoe K, Hasler WL, Lin HC, Maurer AH, et al. Consensus recommendations for gastric emptying scintigraphy: a joint report of the American Neurogastroenterology and Motility Society and the Society of Nuclear Medicine. Am J Gastroenterol. 2008:103 (3):753-63.
- 74. Edwards ST, Cocjin J, Theut SB, Rivard D, Sherman AK, Friesen CA. A comparison of the diagnosis of gastroparesis in 4h pediatric gastric emptying studies versus 2h studies. BMC Gastroenterol. 2019;19 (1):26.
- Chogle A, Saps M. Gastroparesis in children: the benefit of conducting 4-hour scintigraphic gastric-emptying studies. J Pediatr Gastroenterol Nutr. 2013;56 (4):439-42.
- Lin E, Connolly LP, Drubach L, Zurakowski D, DiCanzio J, Mitchell K, et al. Effect of early emptying on quantitation and interpretation of liquid gastric emptying studies of infants and young children. J Nucl Med. 2000;41 (4):596–9.
- 77. Ghoos YF, Maes BD, Geypens BJ, Mys G, Hiele MI, Rutgeerts PJ, et al. Measurement of gastric emptying rate of solids by means of a carbon-labeled octanoic acid breath test. Gastroenterology. 1993;104 (6):1640–7.
- Hauser B, Roelants M, De Schepper J, Veereman G, Caveliers V, Devreker T, et al. Gastric emptying of solids in children: reference values for the (13) C-octanoic acid breath test. Neuro Gastroenterol Motil. 2016;28 (10):1480-7.
- Lee JS, Camilleri M, Zinsmeister AR, Burton DD, Choi MG, Nair KS, et al. Toward office-based measurement of gastric emptying in symptomatic diabetics using [¹³C]octanoic acid breath test. Am J Gastroenterol. 2000;95 (10):2751–61.
- Braden B, Adams S, Duan LP, Orth KH, Maul FD, Lembcke B, et al. The [¹³C]acetate breath test accurately reflects gastric emptying of liquids in both liquid and semisolid test meals. Gastroenterology. 1995;108 (4):1048–55.
- Gatti C, di Abriola FF, Dall'Oglio L, Villa M, Franchini F, Amarri S. Is the ¹³C-acetate breath test a valid procedure to analyse gastric emptying in children? J Pediatr Surg. 2000;35 (1):62–5.
- Delbende B, Perri F, Couturier O, Leodolter A, Mauger P, Bridgi B, et al. ¹³C-octanoic acid breath test for gastric emptying measurement. Eur J Gastroenterol Hepatol. 2000;12 (1):85–91.

 Dickman R, Steinmetz A, Bernnstine H, Groshar D, Niv Y. A novel continuous breath test versus scintigraphy for gastric emptying rate measurement. J Clin Gastroenterol. 2011;45 (1):22-5.

- Zahn A, Langhans CD, Hoffner S, Haberkorn U, Rating D, Haass M, et al. Measurement of gastric emptying by ¹³C-octanoic acid breath test versus scintigraphy in diabetics. Z Gastroenterol. 2003;41 (5):383–90.
- 85. Rao SS, Camilleri M, Hasler WL, Maurer AH, Parkman HP, Saad R, et al. Evaluation of gastrointestinal transit in clinical practice: position paper of the American and European Neurogastroenterology and Motility Societies. Neuro Gastroenterol Motil. 2011;23 (1):8–23.
- Odunsi ST, Camilleri M, Szarka LA, Zinsmeister AR. Optimizing analysis of stable isotope breath tests to estimate gastric emptying of solids. Neuro Gastroenterol Motil. 2009;21 (7):706-e38.
- Chew CG, Bartholomeusz FD, Bellon M, Chatterton BE. Simultaneous ¹³C/14C dual isotope breath test measurement of gastric emptying of solid and liquid in normal subjects and patients: comparison with scintigraphy. Nucl Med Rev Cent East Eur. 2003;6 (1):29–33.
- Hauser B, De Schepper J, Caveliers V, Salvatore S, Salvatoni A, Vandenplas Y. Variability of the ¹³C-octanoic acid breath test for gastric emptying of solids in healthy children. Aliment Pharmacol Ther. 2006;23 (9):1315–9.
- Van Den Driessche M. Study of gastrointestinal motility in infants and children using ¹³C breath tests. Leuven University Press; 2001.
- Szarka LA, Camilleri M, Vella A, Burton D, Baxter K, Simonson J, et al. A stable isotope breath test with a standard meal for abnormal gastric emptying of solids in the clinic and in research. Clin Gastroenterol Hepatol. 2008;6 (6):635-43.
- Braden B, Peterknecht A, Piepho T, Schneider A, Caspary WF, Hamscho N, et al. Measuring gastric emptying of semisolids in children using the ¹³C-acetate breath test: a validation study. Dig Liver Dis. 2004;36 (4):260–4.
- Barbosa L, Vera H, Moran S, Del Prado M, Lopez-Alarcon M. Reproducibility and reliability of the ¹³C-acetate breath test to measure gastric emptying of liquid meal in infants. Nutrition. 2005;21 (3):289–94.
- Sanaka M, Urita Y, Sugimoto M, Yamamoto T, Kuyama Y. Comparison between gastric scintigraphy and the [¹³C]-acetate breath test with Wagner-Nelson analysis in humans. Clin Exp Pharmacol Physiol. 2006;33 (12):1239–43.
- Hauser B, De Schepper J, Caveliers V, Salvatore S, Salvatoni A, Vandenplas Y. Variability of the ¹³C-acetate breath test for gastric emptying of liquids in healthy children. J Pediatr Gastroenterol Nutr. 2006;42 (4):392–7.
- Camilleri M. Clinical practice. Diabetic gastroparesis. N Engl J Med. 2007;356 (8):820–9.
- Ziessman HA, Chander A, Clarke JO, Ramos A, Wahl RL. The added diagnostic value of liquid gastric emptying compared with solid emptying alone. J Nucl Med. 2009;50 (5):726–31.
- 97. van Beek AP, Emous M, Laville M, Tack J. Dumping syndrome after esophageal, gastric or bariatric surgery: pathophysiology, diagnosis, and management. Obes Rev. 2017;18 (1):68–85.
- 98. Bharucha AE, Camilleri M, Veil E, Burton D, Zinsmeister AR. Comprehensive assessment of gastric emptying with a stable isotope breath test. Neuro Gastroenterol Motil. 2013;25 (1):e60-9.
- Keller J, Binnewies U, Rosch M, Juul Holst J, Beglinger C, Andresen V, et al. Gastric emptying and disease activity in inflammatory bowel disease. Eur J Clin Invest. 2015;45 (12):1234–42.
- Bertram F, Andresen V, Layer P, Keller J. Simultaneous noninvasive measurement of liquid gastric emptying and small bowel transit by combined ¹³C-acetate and H2-lactulose breath test. J Breath Res. 2014;8 (4):046007.

- 101. Van Den Driessche M, Peeters K, Marien P, Ghoos Y, Devlieger H, Veereman-Wauters G. Gastric emptying in formula-fed and breastfed infants measured with the ¹³C-octanoic acid breath test. J Pediatr Gastroenterol Nutr. 1999;29 (1):46–51.
- Hauser B, Roelants M, De Schepper J, Veereman G, Caveliers V,
 Devreker T, et al. Gastric emptying of liquids in children. J Pediatr Gastroenterol Nutr. 2016;62 (3):403–8.
- Eradi B, Wright J, Gibbons NJ, Blackshaw PE, Perkins AC, Wakefield J, et al. Validity of ¹³C octanoic acid breath test for measurement of solid meal gastric emptying time in children. J Pediatr Surg. 2006;41 (12):2062-5.
- Keller J, Beglinger C, Holst JJ, Andresen V, Layer P. Mechanisms of gastric emptying disturbances in chronic and acute inflammation of the distal gastrointestinal tract. Am J Physiol Gastrointest Liver Physiol. 2009;297 (5):G861–8.
- 105. Gonlachanvit S, Chey WD, Goodman KJ, Parkman HP. Effect of meal size and test duration on gastric emptying and gastric myoelectrical activity as determined with simultaneous [¹³C]octanoate breath test and electrogastrography in normal subjects using a muffin meal. Dig Dis Sci. 2001;46 (12):2643–50.
- 106. Sanaka M, Yamamoto T, Ishii T, Kuyama Y. The Wagner-Nelson method can generate an accurate gastric emptying flow curve from CO2 data obtained by a ¹³C-labeled substrate breath test. Digestion. 2004;69 (2):71–8.
- 107. Verbeke K. Will the ¹³C-octanoic acid breath test ever replace scintigraphy as the gold standard to assess gastric emptying? Neuro Gastroenterol Motil. 2009;21 (10):1013-6.
- Sanaka M, Yamamoto T, Osaki Y, Kuyama Y. Assessment of the gastric emptying velocity by the ¹³C-octanoate breath test: deconvolution versus a Wagner-Nelson analysis. J Gastroenterol. 2006;41 (7):638–46.
- Maes BD, Ghoos YF, Geypens BJ, Hiele MI, Rutgeerts PJ. Relation between gastric emptying rate and energy intake in children compared with adults. Gut. 1995;36 (2):183–8.
- Maes BD, Ghoos YF, Geypens BJ, Hiele MI, Rutgeerts PJ. Influence of octreotide on the gastric emptying of solids and liquids in normal healthy subjects. Aliment Pharmacol Ther. 1995;9 (1):11–8.
- 111. Perri F, Pastore M, Zicolella A, Annese V, Quitadamo M, Andriulli A. Gastric emptying of solids is delayed in celiac disease and normalizes after gluten withdrawal. Acta Paediatrica (Oslo). 2000;89 (8):921-5.
- 112. Staelens S, Van den Driessche M, Barclay D, Carrié-Faessler AL, Haschke F, Verbeke K, et al. Gastric emptying in healthy newborns fed an intact protein formula, a partially and an extensively hydrolysed formula. Clin Nutr. 2008;27 (2):264–8.
- Machado RS, Yamamoto E, da Silva Patrício FR, Reber M, Kawakami E. Gastric emptying evaluation in children with erosive gastroesophageal reflux disease. Pediatr Surg Int. 2010;26 (5):473-8.
- 114. Perano SJ, Rayner CK, Kritas S, Horowitz M, Donaghue K, Mpundu-Kaambwa C, et al. Gastric emptying is more rapid in adolescents with type 1 diabetes and impacts on postprandial glycemia. J Clin Endocrinol Metab. 2015;100 (6):2248-53.
- 115. Lipp RW, Schnedl WJ, Hammer HF, Kotanko P, Leb G, Krejs GJ. Effects of postprandial walking on delayed gastric emptying and intragastric meal distribution in longstanding diabetics. Am J Gastroenterol. 2000;95 (2):419–24.
- 116. Ikeda T, Inamori M, Fujisawa N, Iwasaki T, Akiyama T, Akimoto K, et al. Effects of body positions on gastric emptying with enteral nutrition: a crossover study using a continuous real time ¹³C breath test (BreathID system). Hepato-gastroenterology. 2008;55 (86-87):1905-7.
- 117. Akimoto K, Inamori M, Iida H, Endo H, Akiyama T, Ikeda T, et al. Does postprandial coffee intake enhance gastric emptying?: a crossover study using continuous real time ¹³C breath test

- (BreathID system). Hepato-gastroenterology. 2009;56 (91-92):918-20.
- Anjiki H, Sanaka M, Kuyama Y. Dual effects of rabeprazole on solidphase gastric emptying assessed by the ¹³C-octanoate breath test. Digestion. 2005;72 (2-3):189-94.
- 119. DiBaise JK, Park FL, Lyden E, Brand RE, Brand RM. Effects of low doses of erythromycin on the ¹³C Spirulina platensis gastric emptying breath test and electrogastrogram: a controlled study in healthy volunteers. Am J Gastroenterol. 2001;96 (7):2041–50.
- 120. Kulik W, van Weissenbruch MM, Menelik N, Cranendonk A, Kneepkens CM, Lafeber HN. Improved use of the [13C]octanoic acid breath test as intra-individual parameter to study the effect of a prokinetic drug on gastric emptying in preterm infants with oral feeding intolerance. J Chromatogr B Biomed Sci Appl. 2001;750 (1):147–53.
- 121. Hoffmeister A, Mayerle J, Beglinger C, Büchler MW, Bufler P, Dathe K, et al. English language version of the S3-consensus guidelines on chronic pancreatitis: definition, aetiology, diagnostic examinations, medical, endoscopic and surgical management of chronic pancreatitis. Z Gastroenterol. 2015;53 (12): 1447–95.
- 122. Thomas PD, Forbes A, Green J, Howdle P, Long R, Playford R, et al. Guidelines for the investigation of chronic diarrhoea, 2nd edition. Gut. 2003;52 (Suppl 5):1–15.
- American Gastroenterological Association medical position statement: guidelines for the evaluation and management of chronic diarrhea. Gastroenterology. 1999;116(6):1461-3.
- 124. Abu-El-Haija M, Uc A, Werlin SL, Freeman AJ, Georgieva M, Jojkić-Pavkov D, et al. Nutritional considerations in pediatric pancreatitis: a position paper from the NASPGHAN Pancreas Committee and ESPGHAN Cystic Fibrosis/Pancreas Working Group. J Pediatr Gastroenterol Nutr. 2018;67 (1):131–43.
- 125. Stern MEH, Palm B, Posselt H-G, Smaczny C. Leitlinien der Gesellschaft für Pädiatrische Gastroenterologie und Ernährung (GPGE): Mukoviszidose (Cystische Fibrose): Ernährung und exokrine Pankreasinsuffizienz. AWMF; 2011. https://www.muko. info/fileadmin/user_upload/angebote/qualitaetsmanagement/LL_ S1_mukoviszidose_ernaehrung_exokrine_pankreasinsuffizienz.pdf.
- 126. Turck D, Braegger CP, Colombo C, Declercq D, Morton A, Pancheva R, et al. ESPEN-ESPGHAN-ECFS guidelines on nutrition care for infants, children, and adults with cystic fibrosis. Clin Nutr. 2016;35 (3):557-77.
- 127. Vantrappen GR, Rutgeerts PJ, Ghoos YF, Hiele MI. Mixed triglyceride breath test: a noninvasive test of pancreatic lipase activity in the duodenum. Gastroenterology. 1989;96 (4):1126–34.
- Perri F, Andriulli A. Mixed" triglyceride breath test: methodological problems and clinical applications. Rev Med Univ Navarra. 1998;42 (2):99–103.
- 129. Sun DY, Jiang YB, Rong L, Jin SJ, Xie WZ. Clinical application of ¹³C-Hiolein breath test in assessing pancreatic exocrine insufficiency. Hepatobiliary & pancreatic diseases international. HBPD Int. 2003;2 (3):449–52.
- 130. Ritz MA, Fraser RJ, Di Matteo AC, Greville H, Butler R, Cmielewski P, et al. Evaluation of the ¹³C-triolein breath test for fat malabsorption in adult patients with cystic fibrosis. J Gastroenterol Hepatol. 2004;19 (4):448–53.
- Löser C, Brauer C, Aygen S, Hennemann O, Fölsch UR. Comparative clinical evaluation of the ¹³C-mixed triglyceride breath test as an indirect pancreatic function test. Scand J Gastroenterol. 1998;33 (3):327–34.
- Domínguez-Muñoz JE, Nieto L, Vilariño M, Lourido MV, Iglesias-García J. Development and diagnostic accuracy of a breath test for pancreatic exocrine insufficiency in chronic pancreatitis. Pancreas. 2016;45 (2):241-7.

- Keller J, Layer P, Bruckel S, Jahr C, Rosien U. ¹³C-mixed triglyceride breath test for evaluation of pancreatic exocrine function in diabetes mellitus. Pancreas. 2014;43 (6):842–8.
- 134. Evenepoel P, Hiele M, Geypens B, Geboes KP, Rutgeerts P, Ghoos Y. ¹³C-egg white breath test: a non-invasive test of pancreatic trypsin activity in the small intestine. Gut. 2000;46 (1):52-7.
- 135. Ishii Y, Kohno T, Ito A, Suzuki S, Kohno T, Takayama T, et al. Evaluation of pancreatic exocrine secretion using ¹³C-dipeptide (benzoyl-L-tyrosyl-[1-(13)C]alanine) breath test: focusing on pancreatoduodenectomy cases. Pancreas. 2007;35 (4):313-9.
- Löser C, Möllgaard A, Aygen S, Hennemann O, Fölsch UR. ¹³C-starch breath test--comparative clinical evaluation of an indirect pancreatic function test. Z Gastroenterol. 1997;35 (3):187-94.
- Siegmund E, Löhr JM, Schuff-Werner P. [The diagnostic validity of non-invasive pancreatic function tests--a meta-analysis]. Z Gastroenterol. 2004;42 (10):1117–28.
- 138. Wejnarska K, Kolodziejczyk E, Ryzko J, Oracz G. Comparison of 72-hour fecal fat quantification and the ¹³C-mixed triglyceride breath test in assessing pancreatic exocrine sufficiency in children with chronic pancreatitis. Dev Period Med. 2016;20 (3):222-7.
- Nakamura H, Murakami Y, Morifuji M, Uemura K, Hayashidani Y, Sudo T, et al. Analysis of fat digestive and absorptive function after subtotal gastrectomy by a ¹³C-labeled mixed triglyceride breath test. Digestion. 2009;80 (2):98–103.
- 140. Nakamura H, Morifuji M, Murakami Y, Uemura K, Ohge H, Hayashidani Y, et al. Usefulness of a ¹³C-labeled mixed triglyceride breath test for assessing pancreatic exocrine function after pancreatic surgery. Surgery. 2009;145 (2):168–75.
- 141. González-Sánchez V, Amrani R, González V, Trigo C, Picó A, de-Madaria E. Diagnosis of exocrine pancreatic insufficiency in chronic pancreatitis: (13)C-mixed triglyceride breath test versus fecal elastase. Pancreatology. 2017;17 (4):580-5.
- 142. Kent DS, Remer T, Blumenthal C, Hunt S, Simonds S, Egert S, et al. ¹³C-Mixed triglyceride breath test and fecal elastase as an indirect pancreatic function test in cystic fibrosis infants. J Pediatr Gastroenterol Nutr. 2018;66 (5):811–5.
- 143. Suzuki M, Tanaka K, Ohtani K, Kitamura K, Kudo T, Shoji H, et al. Estimation of postoperative fat absorption using the ¹³C mixed-triglyceride breath test in children with choledochal cyst. Eur J Pediatr. 2009;168 (1):35–8.
- 144. Lisowska A, Pogorzelski A, Oracz G, Siuda K, Skorupa W, Rachel M, et al. Oral antibiotic therapy improves fat absorption in cystic fibrosis patients with small intestine bacterial overgrowth. J Cyst Fibros. 2011;10 (6):418–21.
- 145. Herzog DC, Delvin EE, Albert C, Marcotte JE, Pelletier VA, Seidman EG. ¹³C-labeled mixed triglyceride breath test (¹³C MTG-BT) in healthy children and children with cystic fibrosis (CF) under pancreatic enzyme replacement therapy (PERT): a pilot study. Clin Biochem. 2008;41 (18):1489-92.
- 146. van Dijk-van Aalst K, Van Den Driessche M, van Der Schoor S, Schiffelers S, van't Westeinde T, Ghoos Y, et al. ¹³C mixed triglyceride breath test: a noninvasive method to assess lipase activity in children. J Pediatr Gastroenterol Nutr. 2001;32 (5):579-85.
- 147. Ling SC, Amarri S, Slater C, Hollman AS, Preston T, Weaver LT. Liver disease does not affect lipolysis as measured with the ¹³C-mixed triacylglycerol breath test in children with cystic fibrosis. J Pediatr Gastroenterol Nutr. 2000;30 (4):368–72.
- 148. Slater C, Preston T, Weaver LT. Improving the specificity of the [¹³C]mixed triacylglycerol breath test by estimating carbon dioxide production from heart rate. Eur J Clin Nutr. 2006;60 (11):1245-52.

 Domínguez-Muñoz JE, Malfertheiner P. Optimized serum pancreolauryl test for differentiating patients with and without chronic pancreatitis. Clin Chem. 1998;44 (4):869–75.

- Hoekstra LT, de Graaf W, Nibourg GA, Heger M, Bennink RJ, Stieger B, et al. Physiological and biochemical basis of clinical liver function tests: a review. Ann Surg. 2013;257 (1):27–36.
- Dietrich CG, Götze O, Geier A. Molecular changes in hepatic metabolism and transport in cirrhosis and their functional importance. World J Gastroenterol. 2016;22 (1):72–88.
- Afolabi P, Wright M, Wootton SA, Jackson AA. ¹³C-aminopyrine demethylation is decreased in cirrhotic patients with normal biochemical markers. Isot Environ Health Stud. 2013;49 (3):346–56.
- 153. Saadeh S, Behrens PW, Parsi MA, Carey WD, Connor JT, Grealis M, et al. The utility of the ¹³C-galactose breath test as a measure of liver function. Aliment Pharmacol Ther. 2003;18 (10):995–1002.
- 154. Giannini EG, Fasoli A, Borro P, Botta F, Malfatti F, Fumagalli A, et al. ¹³C-galactose breath test and ¹³C-aminopyrine breath test for the study of liver function in chronic liver disease. Clin Gastroenterol Hepatol. 2005;3 (3):279-85.
- 155. Braden B, Faust D, Sarrazin U, Zeuzem S, Dietrich CF, Caspary WF, et al. ¹³C-methacetin breath test as liver function test in patients with chronic hepatitis C virus infection. Aliment Pharmacol Ther. 2005;21 (2):179–85.
- 156. Banasch M, Ellrichmann M, Tannapfel A, Schmidt WE, Goetze O. The non-invasive (13)C-methionine breath test detects hepatic mitochondrial dysfunction as a marker of disease activity in non-alcoholic steatohepatitis. Eur J Med Res. 2011;16 (6):258-64.
- Portincasa P, Grattagliano I, Lauterburg BH, Palmieri VO, Palasciano G, Stellaard F. Liver breath tests non-invasively predict higher stages of non-alcoholic steatohepatitis. Clin Sci (Lond). 2006;111 (2):135–43.
- 158. Rocco A, de Nucci G, Valente G, Compare D, D'Arienzo A, Cimino L, et al. ¹³C-aminopyrine breath test accurately predicts long-term outcome of chronic hepatitis C. J Hepatol. 2012;56 (4):782–7.
- 159. Tribonias G, Margariti E, Tiniakos D, Pectasides D, Papatheodoridis GV. Liver function breath tests for differentiation of steatohepatitis from simple fatty liver in patients with nonalcoholic Fatty liver disease. J Clin Gastroenterol. 2014;48 (1):59–65.
- 160. Korkmaz H, Unler GK, Gokturk HS, Schmidt WE, Kebapcilar L. Noninvasive estimation of disease activity and liver fibrosis in nonalcoholic fatty liver disease using anthropometric and biochemical characteristics, including insulin, insulin resistance, and ¹³C-methionine breath test. Eur J Gastroenterol Hepatol. 2015;27 (10):1137–43.
- 161. Afolabi PR, Scorletti E, Smith DE, Almehmadi AA, Calder PC, Byrne CD. The characterisation of hepatic mitochondrial function in patients with non-alcoholic fatty liver disease (NAFLD) using the (13) C-ketoisocaproate breath test. J Breath Res. 2018;12 (4):046002.
- 162. Banasch M, Emminghaus R, Ellrichmann M, Schmidt WE, Goetze O. Longitudinal effects of hepatitis C virus treatment on hepatic mitochondrial dysfunction assessed by C-methionine breath test. Aliment Pharmacol Ther. 2008;28 (4):443-9.
- 163. Fierbinteanu-Braticevici C, Plesca DA, Tribus L, Panaitescu E, Braticevici B. The role of ¹³C-methacetin breath test for the non-invasive evaluation of nonalcoholic fatty liver disease. J Gastro-intestin Liver Dis. 2013;22 (2):149–56.
- 164. Kempinski R, Neubauer K, Wieczorek S, Dudkowiak R, Jasinska M, Poniewierka E. ¹³C-Methacetin breath testing in patients with nonalcoholic fatty liver disease. Adv Clin Exp Med. 2016;25 (1):77–81.
- 165. Dinesen L, Caspary WF, Chapman RW, Dietrich CF, Sarrazin C, Braden B. ¹³C-methacetin-breath test compared to also noninvasive biochemical blood tests in predicting hepatic fibrosis and cirrhosis in chronic hepatitis C. Digestive and liver disease. Dig Liver Dis. 2008;40 (9):743–8.

166. Giannini EG, Savarino V. Relationship between ¹³C-aminopyrine breath test and the MELD score and its long-term prognostic use in patients with cirrhosis. Dig Dis Sci. 2013;58 (10):3024–8.

- 167. Park GJ, Katelaris PH, Jones DB, Seow F, Le Couteur DG, Ngu MC. Validity of the ¹³C-caffeine breath test as a noninvasive, quantitative test of liver function. Hepatology. 2003;38 (5):1227–36.
- 168. Ishii T, Furube M, Hirano S, Takatori K, Iida K, Kajiwara M. Evaluation of ¹³C-phenylalanine and ¹³C-tyrosine breath tests for the measurement of hepatocyte functional capacity in patients with liver cirrhosis. Chem Pharm Bull. 2001;49 (12):1507–11.
- Lykke Eriksen P, Sørensen M, Grønbæk H, Hamilton-Dutoit S, Vilstrup H, Thomsen KL. Non-alcoholic fatty liver disease causes dissociated changes in metabolic liver functions. Clin Res Hepatol Gastroenterol. 2019;43 (5):551–60.
- 170. Banasch M, Goetze O, Hollborn I, Hochdorfer B, Bulut K, Schlottmann R, et al. ¹³C-methionine breath test detects distinct hepatic mitochondrial dysfunction in HIV-infected patients with normal serum lactate. J Acquir Immune Defic Syndr. 2005;40 (2):149–54.
- 171. Forestier J, Dumortier J, Guillaud O, Ecochard M, Roman S, Boillot O, et al. Noninvasive diagnosis and prognosis of liver cirrhosis: a comparison of biological scores, elastometry, and metabolic liver function tests. Eur J Gastroenterol Hepatol. 2010;22 (5):532–40.
- Palmieri VO, Grattagliano I, Minerva F, Pollice S, Palasciano G,
 Portincasa P. Liver function as assessed by breath tests in patients
 with hepatocellular carcinoma. J Surg Res. 2009;157 (2):199–207.
- 173. Goetze O, Breuer M, Geier A, Fried M, Weber A, Jochum W, et al. The ¹³C-methactin breath test is non-inferior to liver biopsy in predicting liver-related death and transplantation: a 7-year prospective follow-up study in 132 patients with chronic hepatitis C infection. GastroHep. 2020;2 (6):344–50.
- 174. Barzakova ES, Schulze-Hagen M, Zimmermann M, Lurje G, Bednarsch J, Pedersoli F, et al. Monitoring liver function of patients undergoing transarterial chemoembolization (TACE) by a ¹³C breath test (LiMAx). Cardiovasc Intervent Radiol. 2019;42 (12):1702–8.
- 175. Makridis G, Oldhafer KJ. First intraoperative measurement of liver functional capacity during liver surgery using the (13) Cmethacetin breath test: early results of a pilot study. J Hepato-Biliary-Pancreatic Sci. 2020;27 (5):280-1.
- 176. Jara M, Reese T, Malinowski M, Valle E, Seehofer D, Puhl G, et al. Reductions in post-hepatectomy liver failure and related mortality after implementation of the LiMAx algorithm in preoperative work-up: a single-centre analysis of 1170 hepatectomies of one or more segments. HPB (Oxford). 2015;17 (7):651–8.
- Stockmann M, Lock JF, Malinowski M, Niehues SM, Seehofer D, Neuhaus P. The LiMAx test: a new liver function test for predicting postoperative outcome in liver surgery. HPB (Oxford). 2010;12 (2):139-46.
- 178. Duro D, Fitzgibbons S, Valim C, Yang CF, Zurakowski D, Dolan M, et al. [13C]Methionine breath test to assess intestinal failure-associated liver disease. Pediatr Res. 2010;68 (4):349–54.
- 179. Shteyer E, Lalazar G, Hemed N, Pappo O, Granot E, Yerushalmi B, et al. Continuous ¹³C-methacetin breath test differentiates biliary atresia from other causes of neonatal cholestasis. J Pediatr Gastroenterol Nutr. 2013;56 (1):60–5.
- Wada M, Wada Y, Uchiyama M, Kajiwara M, Takatori K. (13)Cphenylalanine breath test correlates with liver fibrosis in postoperative biliary atresia. Pediatr Int. 2007;49 (6):836-41.
- Turki A, Murthy G, Ueda K, Cheng B, Giezen A, Stockler-Ipsiroglu S, et al. Minimally invasive (13)C-breath test to examine phenylalanine metabolism in children with phenylketonuria. Mol Genet Metabol. 2015;115 (2-3):78–83.
- Lynch T, Price A. The effect of cytochrome P450 metabolism on drug response, interactions, and adverse effects. Am Fam Physician. 2007;76 (3):391–6.

- 183. Herold C, Ganslmayer M, Ocker M, Zopf S, Gailer B, Hahn EG, et al. Inducibility of microsomal liver function may differentiate cirrhotic patients with maintained compared with severely compromised liver reserve. J Gastroenterol Hepatol. 2003;18 (4):445-9.
- 184. Jonderko K, Skałba P, Kasicka-Jonderko A, Kamińska M, Bizior-Frymus D, Dyja R. Impact of combined oral contraceptives containing ethinylestradiol on the liver microsomal metabolism. Eur J Contracept Reprod Health Care. 2013;18 (4):284–92.
- 185. Banerjee D, Vikram N, Mishra P, Bhatt R, Prakash S, Misra A. Correlation of a [13C]glucose breath test with surrogate markers of insulin resistance in urban and rural Asian Indians. Metab Syndr Relat Disord. 2009;7 (3):215-9.
- 186. Ghosh C, Maity A, Banik GD, Som S, Chakraborty A, Selvan C, et al. Non-invasive ¹³C-glucose breath test using residual gas analyzer-mass spectrometry: a novel tool for screening individuals with pre-diabetes and type 2 diabetes. J Breath Res. 2014;8 (3): 036001.
- 187. Hussain M, Jangorbhani M, Schuette S, Considine RV, Chisholm RL, Mather KJ. [13C]glucose breath testing provides a noninvasive measure of insulin resistance: calibration analyses against clamp studies. Diabetes Technol Therapeut. 2014;16 (2):102–12.
- Jetha MM, Nzekwu U, Lewanczuk RZ, Ball GD. A novel, noninvasive ¹³C-glucose breath test to estimate insulin resistance in obese prepubertal children. J Pediatr Endocrinol Metab. 2009;22 (11):1051-9.
- 189. Kawagoe N, Kano O, Kijima S, Tanaka H, Takayanagi M, Urita Y. Investigation of metabolism of exogenous glucose at the early stage and onset of diabetes mellitus in Otsuka long-evans Tokushima fatty rats using [1, 2, 3-13C]glucose breath tests. PloS One. 2016;11 (8):e0160177.
- Lewanczuk RZ, Paty BW, Toth EL. Comparison of the [¹³C] glucose breath test to the hyperinsulinemic-euglycemic clamp when determining insulin resistance. Diabetes Care. 2004;27 (2): 441–7
- 191. Maldonado-Hernandez J, Martinez-Basila A, Salas-Fernandez A, Navarro-Betancourt JR, Pina-Aguero MI, Bernabe-Garcia M. The ¹³C-glucose breath test for insulin resistance assessment in adolescents: comparison with fasting and post-glucose stimulus surrogate markers of insulin resistance. J Clin Res Pediatr Endocrinol. 2016;8 (4):419-24.
- 192. Mizrahi M, Lalazar G, Adar T, Raz I, Ilan Y. Assessment of insulin resistance by a ¹³C glucose breath test: a new tool for early diagnosis and follow-up of high-risk patients. Nutr J. 2010:9:25.
- 193. Salas-Fernandez A, Maldonado-Hernandez J, Martinez-Basila A, Martinez-Razo G, Jasso-Saavedra F. The ¹³C-glucose breath test is a valid non-invasive screening tool to identify metabolic syndrome in adolescents. Clin Chem Lab Med. 2015;53 (1):133–8.
- 194. Takemoto I, Kawagoe N, Kijima S, Sasaki Y, Watanabe T, Urita Y. (13)C-glucose breath tests: a non-invasive method for detecting early clinical manifestations of exogenous glucose metabolism in type 2 diabetic patients. Acta Diabetol. 2019;56 (4):449-56.
- 195. Berthold HK, Schober P, Scheurlen C, Marklein G, Horre R, Gouni-Berthold I, et al. Use of the lactose-[¹³C]ureide breath test for diagnosis of small bowel bacterial overgrowth: comparison to the glucose hydrogen breath test. J Gastroenterol. 2009;44 (9):944–51.
- Wutzke KD, Glasenapp B. The use of ¹³C-labelled glycosyl ureides for evaluation of orocaecal transit time. Eur J Clin Nutr. 2004;58 (4):568–72.
- Morrison DJ, Dodson B, Preston T, Weaver LT. Gastrointestinal handling of glycosyl [¹³C]ureides. Eur J Clin Nutr. 2003;57 (8):1017-24.

- 198. Van Den Driessche M, Van Malderen N, Geypens B, Ghoos Y, Veereman-Wauters G. Lactose-[13C]ureide breath test: a new, noninvasive technique to determine orocecal transit time in children. J Pediatr Gastroenterol Nutr. 2000;31 (4):433–8.
- 199. Jonderko K, Kasicka-Jonderko A, Syrkiewicz-Trepiak D, Blonska-Fajfrowska B. Feasibility of a breath test with a substrate of natural ¹³C-abundance and isotope-selective non-dispersive infrared spectrometry: a preliminary study. J Gastroenterol Hepatol. 2005;20 (8):1228-34.
- Robayo-Torres CC, Diaz-Sotomayor M, Hamaker BR, Baker SS, Chumpitazi BP, Opekun AR, et al. ¹³C-Labeled-Starch breath test in congenital sucrase-isomaltase deficiency. J Pediatr Gastroenterol Nutr. 2018;66 (Suppl 3):S61.
- Robayo-Torres CC, Opekun AR, Quezada-Calvillo R, Villa X, Smith EO, Navarrete M, et al. ¹³C-breath tests for sucrose digestion in congenital sucrase isomaltase-deficient and sacrosidase-supplemented patients. J Pediatr Gastroenterol Nutr. 2009;48 (4):412-8.
- 202. de Wildt SN, Berns MJ, van den Anker JN. ¹³C-erythromycin breath test as a noninvasive measure of CYP3A activity in newborn infants: a pilot study. Ther Drug Monit. 2007;29 (2):225–30.
- Furuta T, Iwaki T, Umemura K. [¹³C]pantoprazole breath test as a predictor of the anti-platelet function of clopidogrel. Eur J Clin Pharmacol. 2010;66 (5):457–63.
- 204. Furuta T, Kodaira C, Nishino M, Yamade M, Sugimoto M, Ikuma M, et al. [¹³C]-pantoprazole breath test to predict CYP2C19 phenotype and efficacy of a proton pump inhibitor, lansoprazole. Aliment Pharmacol Ther. 2009;30 (3):294–300.
- Modak AS, Klyarytska I, Kriviy V, Tsapyak T, Rabotyagova Y. The
 effect of proton pump inhibitors on the CYP2C19 enzyme activity
 evaluated by the pantoprazole-(13)C breath test in GERD patients:
 clinical relevance for personalized medicine. J Breath Res. 2016;10
 (4):046017.
- 206. Leeder JS, Pearce RE, Gaedigk A, Modak A, Rosen DI. Evaluation of a [13C]-dextromethorphan breath test to assess CYP2D6 phenotype. J Clin Pharmacol. 2008;48 (9):1041–51.
- 207. Opdam FL, Dezentje VO, den Hartigh J, Modak AS, Vree R, Batman E, et al. The use of the ¹³C-dextromethorphan breath test for phenotyping CYP2D6 in breast cancer patients using tamoxifen: association with CYP2D6 genotype and serum endoxifen levels. Canc Chemother Pharmacol. 2013;71 (3):593–601.
- 208. Safgren SL, Suman VJ, Kosel ML, Gilbert JA, Buhrow SA, Black JL, et al. Evaluation of CYP2D6 enzyme activity using a ¹³C-dextromethorphan breath test in women receiving adjuvant tamoxifen. Pharmacogenetics Genom. 2015;25 (4):157–63.
- Mattison LK, Ezzeldin H, Carpenter M, Modak A, Johnson MR, Diasio RB. Rapid identification of dihydropyrimidine dehydrogenase deficiency by using a novel 2-¹³C-uracil breath test. Clin Canc Res. 2004;10 (8):2652-8.
- Mattison LK, Fourie J, Hirao Y, Koga T, Desmond RA, King JR, et al.
 The uracil breath test in the assessment of dihydropyrimidine dehydrogenase activity: pharmacokinetic relationship between expired ¹³CO2 and plasma [2-13C]dihydrouracil. Clin Canc Res. 2006;12 (2):549-55.

How to cite this article: Keller J, Hammer HF, Afolabi PR, Benninga M, Borrelli O, Dominguez-Munoz E, et al. European guideline on indications, performance and clinical impact of ¹³C-breath tests in adult and pediatric patients: An EAGEN, ESNM, and ESPGHAN consensus, supported by EPC. United European Gastroenterol J. 2021;1-29. https://doi.org/10.1002/ueg2.12099